



200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (sterile)



Product Highlights

- **Graduated** at 10 µl, 50 µl and 100 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified** RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	S1113-1810
Pack Size	960 Tips (10 Racks × 96 Tips)
Color Code Wafer	Yellow
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Tip Type	UltraPoint®
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with	<p>Abimed: BioPette-E, Discovery Comfort, LABMATE Optima</p> <p>Biohit®: eLine, ePet, mLine, Picus, Proline</p> <p>BRAND®: Transferpette, Transferpette Electronic, Transferpette S</p> <p>Capp®: CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel</p> <p>Eppendorf®: Reference, Reference 2, Research, Research plus, Research pro, Xplorer</p> <p>Gilson®: Pipetman Concept, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p>Hamilton®: SoftGrip Pipette, SofTouch</p> <p>HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p>Nichiryo®: 00-NP7, 50F, 50V, NLE, NPLO, NPX</p> <p>Rainin®: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p>Sartorius®: eLine, ePet, mLine, Picus, Picus Nxt, Proline, Proline Plus, Tacta</p> <p>Socorex®: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital</p> <p>STARLAB®: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p>Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finn Focus Univ., Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Turquoise, Finnpipette Novus Yellow</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (sterile)

Improved by you - New System Components Overview

The perfect rack

- › Simple! The system comprises only two racks (one small, one large) for less complexity
- › Small rack for use with tip volumes 10 µl to 300 µl; large rack for use with tip volumes 1000 µl and 1250µl XL
- › Quick and easy identification of tips by colour-coded wafers with laser printed information and transparent lid
- › Stackable
- › Reusable
- › Autoclavable

The clever refill

- › The new refill is divided into two separate, individually-wrapped stacks of five wafers
- › Completely sealed, the tips are safe from contamination by the tight fitting top shell and unique spacers
- › Economical, space-saving design
- › Autoclavable when racked

The self-standing bag

- › New, self-standing bag for easy and convenient access and storage
- › Tamper-proof, easy-open perforation strip
- › Resealable zip-lock for contamination-free storage
- › Recyclable polyethylene bag
- › Economical – ideal for teaching labs or non-critical work

The compact stack rack


- › Next generation design with wide base for great stability
- › Transparent lid covers the tips when not in use
- › 960 tips for high throughput use, ideal for multichannel users
- › The complete stack rack is autoclavable and recyclable

The unique filter refill













- › A brand new, economical and ecofriendly TipOne® filter tip refill system
- › Fast, easy and safe: The filter tips are refilled with the top and bottom shell intact, thus providing a fast and completely protected refill system that will not compromise sterility
- › Using sterile filter refills to reload racks directly during periods of long use in hoods and PCR-workstations saves valuable space
- › Not only does the refill system generate less waste, the bottom shells are nestable so your waste takes up less storage space too

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk















All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Rack (non-sterile) Sterile: No Pack Size: 960 Tips (10 Rack × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Rack × 96 Tips)	S1113-1800
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1113-1800-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1113-1810
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1113-1810-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1113-1700
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1113-1700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl UltraPoint® graduated TipOne® Tip, Natural Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1113-1710
 200 µl UltraPoint® graduated TipOne® Tip, Natural Refills (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Refills × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Refills × 96 Tips)	S1113-1710-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)  Create PDF Data Sheet	960 Tips (1 Stack rack × 960 Tips)	S1113-1200
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)  Create PDF Data Sheet	9,600 Tips (10 Stack racks × 960 Tips)	S1113-1200-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	S1113-1000
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 25,000 Tips (25 Bags × 1000 Tips)  Create PDF Data Sheet	25,000 Tips (25 Bags × 1000 Tips)	S1113-1000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk