



Racks for 15ml centrifuge tubes (15-place)



Product Highlights

- Chemically-resistant epoxy-coated wire racks
- Suitable for use in the hood
- Temperature use: -80 °C to 121 °C
- Dimensions (WxDxH): 127 x 76 x 176 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2396-1515
Pack Size	1 Piece (1 Rack × 1 Piece)
Dimensions (W × D × H)	127 mm × 76 mm × 176 mm
Dimensions (W × D)	127 mm × 76 mm
External width	127 mm
Max. tube diameter	20 mm
Autoclavable	Yes
Tube Style	15 ml Centrifuge
Color	Blue
Operating temperature	-80°C to 121°C
Max. centrifuge tubes or screwcap microtubes	15ml tubes
Places	15

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Racks for 15ml centrifuge tubes (15-place)













Chemically-resistant epoxy-coated wire racks, suitable for use in the hood. Temperature use: -80 °C to 121 °C. Autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk







All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Racks for 13mm blood collection tubes (72-place) Tube Style: 13 mm Blood Collection Places: 72  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-7213
	Racks for 16mm blood collection tubes (72-place) Tube Style: 16 mm Blood Collection Places: 72  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-7216
	Racks for 30 ml universal containers (50-Place) Tube Style: 30 ml Universal Places: 50  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1750
	Racks for 30 ml universal containers (12-Place) Tube Style: 30 ml Universal Places: 12  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1712
	Racks for 15ml centrifuge tubes (15-place) Tube Style: 15 ml Centrifuge Places: 15  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1515
	Racks for 50 ml centrifuge tubes (8-place) Tube Style: 50 ml Centrifuge Places: 8  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1550

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Racks for 7.0 ml bijou containers (12-Place) Tube Style: 7.0 ml Bijou Places: 12  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1612
 Racks for 7.0 ml bijou containers (50-Place) Tube Style: 7.0 ml Bijou Places: 50  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1650

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk