



9 Litre Oblong Cooler Container, Red (521 x 295 x 165 mm)



Product Highlights

- The snug-fitting lid minimises thawing and evaporation
- Carrying handles
- Handy pour spout for easy emptying
- Tapered sides allow stacking when not in use
- Temperature use down to -80 °C

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2396-6004
Pack Size	1 Piece (1 Box × 1 Piece)
Dimensions (W × D × H)	521 mm × 295 mm × 165 mm
Dimensions (W × D)	521 mm × 295 mm
External width	521 mm
Autoclavable	No
Volume	9 l
Material	Polyurethane
Color	Red

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk









More informations about 9 Litre Oblong Cooler Container, Red (521 x 295 x 165 mm)

Containers for temperature sensitive samples

Rigid moulded polyurethane containers for chilling your temperature sensitive samples. The Cooler Containers provide excellent insulation properties for very low temperature use with ice and ice-salt slurries.



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	4 Litre Oblong Cooler Container, Blue (380 x 250 x 140 mm) Volume: 4 l  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2396-6001
	1 Litre Square Cooler Container, Green (270 x 195 x 115 mm) Volume: 1 l  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2396-6002
	9 Litre Oblong Cooler Container, Red (521 x 295 x 165 mm) Volume: 9 l  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2396-6004

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk