



## 2.0 ml Crystal Clear Microcentrifuge Tube



### Product Highlights

- Ultra smooth tube bottoms for 100% visibility
- Thin cap membrane for easy sample access with a syringe needle
- Max. centrifugation rate: 14,000–20,000 x g
- Temperature range: -80°C to 121°C
- Autoclavable
- Certified RNase, DNase and Pyrogen free

## Product Applications

Sample preparation

Centrifugation

Every application where a penetration of the cap is necessary

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E1420-2000
Pack Size	500 Tubes (1 Bags × 500 Tubes)
Sterile	No
Autoclavable	Yes
Volume	2.0 ml
Material	Polypropylene
DNA free	No
DNase free	Yes
Free of endotoxins	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Crystal clear
g-Force max.	14,000 x g - 20,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 2.0 ml Crystal Clear Microcentrifuge Tube













- › Crystal Clear polypropylene tubes for the ultimate in sample visibility
- › Ultra smooth tube bottoms for 100% visibility
- › Thin cap membrane for easy sample access with a syringe needle
- › Frosted writing surface
- › Max. centrifugation rate: 14,000–20,000 x g
- › Temperature range: -80°C to 121°C.
- › Autoclavable
- › Certified RNase, DNase and Pyrogen free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations


PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.6 ml Crystal Clear Microcentrifuge Tube</a> <b>Volume:</b> 0.6 ml <b>Sterile:</b> No  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (2 Bags × 500 Tubes)	E1405-0600
 <a href="#">0.6 ml Crystal Clear Microcentrifuge Tube</a> <b>Volume:</b> 0.6 ml <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	500 Tubes (5 Bags × 100 Tubes)	E1405-0610
 <a href="#">1.5 ml Crystal Clear Microcentrifuge Tube</a> <b>Volume:</b> 1.5 ml <b>Sterile:</b> No  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bags × 500 Tubes)	E1415-1500
 <a href="#">1.5 ml Crystal Clear Microcentrifuge Tube</a> <b>Volume:</b> 1.5 ml <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-1510
 <a href="#">2.0 ml Crystal Clear Microcentrifuge Tube</a> <b>Volume:</b> 2.0 ml <b>Sterile:</b> No  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bags × 500 Tubes)	E1420-2000
 <a href="#">2.0 ml Crystal Clear Microcentrifuge Tube</a> <b>Volume:</b> 2.0 ml <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1420-2010

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes</a>	100 Pcs. (1 Bag × 100 Pcs.)	I1415-1508

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)