




## 50 ml Centrifuge Tube, Conical, Loose (Sterile)



### Product Highlights

- Clear, pre-sterilized polypropylene tubes with blue HDPE screw caps
- White writing surface
- Clearly printed graduations
- Certified RNase, DNase, DNA and Pyrogen free
- Non-Cytotoxic
- Centrifugation Rate Conical: Swing Bucket Rotor: 3,200 x g (swing bucket) / 9,500 x g (fixed angle rotors)

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E1450-0200
Pack Size	500 Tubes (10 Bags × 50 Tubes)
Dimensions (W × H)	34.9 mm × 115.5 mm
Dimensions (W × L)	34.9 mm
External width	34.9 mm
Diameter	23 mm
Sterile	Yes
Volume	50 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	No
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Clear
g-Force max.	3,200 × g (swing bucket) / 9,500 × g (fixed angle rotors)
Non-cytotoxic	Yes
Tube Base	Conical

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 50 ml Centrifuge Tube, Conical, Loose (Sterile)

Don't get in a spin! Use Starlab centrifuge tubes!

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)