



## 1.8 ml Cryovial Vials with External Thread, Self-Sealing Cap Conical (Sterile)



### Product Highlights

- External Threaded Tube / Self-Sealing Cap
- White writing area and highly visible black graduations at 100 µl increments
- Sterile (SAL 10-□)
- Certified RNase, DNase, DNA, PCR inhibitors and Endotoxin free
- Cap can be colour-coded using the coloured cap inserts


## Product Applications

Storage and transportation of biological material

Ideal for vapor-phase liquid nitrogen freezing (Not suitable for immersion)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                       |   |
|-----------------------|---|
| Art. No.              | E3090-6321  |
| Pack Size             | 500 Pcs. (10 Bags × 50 Pcs.)                              |
| Sterile               | Yes   |
| Volume                | 1.8 ml  |
| Cap style             | External Thread, Self-Sealing Cap                         |
| Material              | Polypropylene   |
| DNA free              | Yes   |
| DNase free            | Yes   |
| PCR inhibitor free    | Yes   |
| RNase free            | Yes   |
| Color                 | Natural   |
| Operating temperature | -196°C (only be used in the gas phase of liquid nitrogen) |
| Tube Base             | Conical   |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 1.8 ml Cryovial Vials with External Thread, Self-Sealing Cap Conical (Sterile)









Starlab's Cryogenic Vials have been developed for the storage and transportation of biological material.

- › Made from 100 % polypropylene
- › For storage down to -196 °C; should only be used in the gas phase of liquid nitrogen
- › White writing area and highly-visible black graduations at 100 µl increments
- › Sterile (SAL 10-6)
- › Certified RNase, DNase, DNA, PCR Inhibitor and Endotoxin free
- › The internal threaded vials with a silicone seal in cap provide the highest possible impermeability
- › Caps can be colour coded with the Coloured Cap Inserts

Note: this cryovial may not fit in all storage canes due to its diameter.



## All Variations

| PRODUCT NAME   | PACKAGING SIZE               | ART. NO.   |
|--|------------------------------|------------|
| <br><a href="#">1.0 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile)</a><br><b>Volume:</b> 1.0 ml<br><b>Tube Base:</b> Skirted<br> <a href="#">Create PDF Data Sheet</a>      | 500 Pcs. (10 Bags × 50 Pcs.) | E3090-6212 |
| <br><a href="#">1.8 ml Cryovial Vials with External Thread, Self Sealing Seal Cap, Skirted (Sterile)</a><br><b>Volume:</b> 1.8 ml<br><b>Tube Base:</b> Skirted<br> <a href="#">Create PDF Data Sheet</a> | 500 Pcs. (10 Bags × 50 Pcs.) | E3090-6222 |
| <br><a href="#">4.5 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile)</a><br><b>Volume:</b> 4.5 ml<br><b>Tube Base:</b> Skirted<br> <a href="#">Create PDF Data Sheet</a>    | 300 Pcs. (6 Bags × 50 Pcs.)  | E3090-6242 |
| <br><a href="#">1.8 ml Cryovial Vials with External Thread, Self-Sealing Cap Conical (Sterile)</a><br><b>Volume:</b> 1.8 ml<br><b>Tube Base:</b> Conical<br> <a href="#">Create PDF Data Sheet</a>   | 500 Pcs. (10 Bags × 50 Pcs.) | E3090-6321 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

| PRODUCT NAME  |   | PACKAGING SIZE                     | ART. NO.   |
|---|---|------------------------------------|------------|
|  | <a href="#">9.5 mm diameter, Cryogenic Spot Labels, Mixed Colours</a> | 5,000 Labels (1 Box × 5000 Labels) | E9189-0099 |
|  | <a href="#">StarTag Cryo (38 x 15 mm), White</a>                      | 1,000 Labels (1 Box × 1000 Labels) | E9189-3009 |
|  | <a href="#">Storage Canes for 5 x 1.2 / 2.0 ml Cryogenic Vials</a>    | 12 Pcs. (1 Pack × 12 Pcs.)         | E2005-5000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)