



1.8 ml Cryovial Vials with External Thread, Self Sealing Seal Cap, Skirted (Sterile)



Product Highlights

- External Threaded Tube / Self-Sealing Cap
- White writing area and highly visible black graduations at 100 µl increments
- Sterile (SAL 10-□)
- Certified RNase, DNase, DNA, PCR inhibitors and Endotoxin free
- Cap can be colour-coded using the coloured cap inserts

Product Applications

Storage and transportation of biological material

Ideal for vapor-phase liquid nitrogen freezing (Not suitable for immersion)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E3090-6222
Pack Size	500 Pcs. (10 Bags × 50 Pcs.)
Sterile	Yes
Volume	1.8 ml
Cap style	External Thread, Self-Sealing Cap
Material	Polypropylene
DNA free	Yes
DNase free	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural
Operating temperature	-196°C (only be used in the gas phase of liquid nitrogen)
Tube Base	Skirted

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 1.8 ml Cryovial Vials with External Thread, Self Sealing Seal Cap, Skirted (Sterile)

Starlab's Cryogenic Vials have been developed for the storage and transportation of biological material.









- › Made from 100 % polypropylene
- › For storage down to -196 °C; should only be used in the gas phase of liquid nitrogen
- › White writing area and highly-visible black graduations at 100 µl increments
- › Sterile (SAL 10-6)
- › Certified RNase, DNase, DNA, PCR Inhibitor and Endotoxin free
- › The internal threaded vials with a silicone seal in cap provide the highest possible impermeability
- › The skirted vials lock into racks for easy single-handed work
- › Caps can be colour coded with the Coloured Cap Inserts

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1.0 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile) Volume: 1.0 ml Tube Base: Skirted  Create PDF Data Sheet	500 Pcs. (10 Bags × 50 Pcs.)	E3090-6212
 1.8 ml Cryovial Vials with External Thread, Self Sealing Seal Cap, Skirted (Sterile) Volume: 1.8 ml Tube Base: Skirted  Create PDF Data Sheet	500 Pcs. (10 Bags × 50 Pcs.)	E3090-6222
 4.5 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile) Volume: 4.5 ml Tube Base: Skirted  Create PDF Data Sheet	300 Pcs. (6 Bags × 50 Pcs.)	E3090-6242
 1.8 ml Cryovial Vials with External Thread, Self-Sealing Cap Conical (Sterile) Volume: 1.8 ml Tube Base: Conical  Create PDF Data Sheet	500 Pcs. (10 Bags × 50 Pcs.)	E3090-6321

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	9.5 mm diameter, Cryogenic Spot Labels, Mixed Colours	5,000 Labels (1 Box × 5000 Labels)	E9189-0099
	StarTag Cryo (38 x 15 mm), White	1,000 Labels (1 Box × 1000 Labels)	E9189-3009
	Storage Canes for 5 x 1.2 / 2.0 ml Cryogenic Vials	12 Pcs. (1 Pack × 12 Pcs.)	E2005-5000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk