



1.5 ml Screw Cap Tube, Ribbed, Conical, Natural



Product Highlights

• Optically clear

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.









General Data

Art. No.	
	E1415-2210
Pack Size 5	500 Tubes (1 Pack × 500 Tubes)
Sterile	No
Autoclavable	Yes
Volume 1	1.5 ml
Tube Style	Ribbed
Material F	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural
g-Force max.	20,000 x g
Tube Base (Conical

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.



MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







More informations about 1.5 ml Screw Cap Tube, Ribbed, Conical, Natural

Mix and match – tubes and caps!

Polypropylene screw cap tubes that are suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave (121 °C).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

starlab

Starlab (UK) LTD





All Variations

PRODUCT NAME	·	PACKAGING SIZE	ART. NO.
	0.5 ml Screw Cap Tube, Ribbed, Skirted, Natural Volume: 0.5 ml Tube Base: Skirted Tube Style: Ribbed Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1405-2120
	O.5 ml Screw Cap Tube, Plain, Conical, Natural Volume: 0.5 ml Tube Base: Conical Tube Style: Plain Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1405-2130
	0.5 ml Screw Cap Tube, Plain, Skirted, Natural Volume: 0.5 ml Tube Base: Skirted Tube Style: Plain Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1405-2140
	1.5 ml Screw Cap Tube, Ribbed, Conical, Natural Volume: 1.5 ml Tube Base: Conical Tube Style: Ribbed Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2210
	1.5 ml Screw Cap Tube, Ribbed, Skirted, Natural Volume: 1.5 ml Tube Base: Skirted Tube Style: Ribbed Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2220

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.



5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







PRODUCT NAME	i e	PACKAGING SIZE	ART. NO.
	1.5 ml Screw Cap Tube, Graduated, Conical, Natural Volume: 1.5 ml Tube Base: Conical Tube Style: Graduated Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2230
	1.5 ml Screw Cap Tube, Plain, Skirted, Natural Volume: 1.5 ml Tube Base: Skirted Tube Style: Plain Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2240
Ī	2.0 ml Screw Cap Tube, Ribbed, Skirted, Natural Volume: 2.0 ml Tube Base: Skirted Tube Style: Ribbed Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1420-2320
	2.0 ml Screw Cap Tube, Graduated, Conical, Natural Volume: 2.0 ml Tube Base: Conical Tube Style: Graduated Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1420-2330
	2.0 ml Screw Cap Tube, Graduated, Skirted, Natural Volume: 2.0 ml Tube Base: Skirted Tube Style: Graduated Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1420-2340

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.



5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

