



2.0 ml Screw Cap Tube, Graduated, Skirted, Amber



Product Highlights

- Graduated versions marked at every 100 μ l to provide easy sample estimation
- Routinely tested to 20,000 x g
- All tubes and caps are autoclavable
- Certified RNase, DNase and DNA free
- 100% medical grade polypropylene
- Suitable for light sensitive samples

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E1420-2347
Pack Size	500 Tubes (1 Bag × 500 Tubes)
Sterile	No
Autoclavable	Yes
Volume	2.0 ml
Tube Style	Graduated
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Amber
g-Force max.	20,000 × g
Tube Base	Skirted

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk













More informations about 2.0 ml Screw Cap Tube, Graduated, Skirted, Amber

Mix and match – tubes and caps!

Polypropylene screw cap tubes made from 100% medical grade polypropylene. Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave (121 °C). Available with or without graduations, the graduated versions are marked at every 100 µl to provide easy sample estimation. Offered with either skirted or conical bases. Routinely tested to 20,000 x g. Caps are available separately; all tubes and caps are autoclavable. Certified RNase, DNase and DNA free.







All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	0.5 ml Screw Cap Tube, Ribbed, Skirted, Amber Volume: 0.5 ml Tube Style: Ribbed Tube Base: Skirted  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1405-2127
	1.5 ml Screw Cap Tube, Ribbed, Skirted, Amber Volume: 1.5 ml Tube Style: Ribbed Tube Base: Skirted  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1415-2227
	1.5 ml Screw Cap Tube, Graduated, Conical, Amber Volume: 1.5 ml Tube Style: Graduated Tube Base: Conical  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1415-2237
	2.0 ml Screw Cap Tube, Graduated, Conical, Amber Volume: 2.0 ml Tube Style: Graduated Tube Base: Conical  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1420-2337
	2.0 ml Screw Cap Tube, Graduated, Skirted, Amber Volume: 2.0 ml Tube Style: Graduated Tube Base: Skirted  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1420-2347

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1.5 ml Screw Cap Tube, Plain, Skirted, Amber Volume: 1.5 ml Tube Style: Plain Tube Base: Skirted  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1415-2247
 0.5 ml Screw Cap Tube, Plain, Skirted, Amber Volume: 0.5 ml Tube Style: Plain Tube Base: Skirted  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	E1405-2147

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk