




## 1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile), Case



### Product Highlights

- Graduated at 100 µl, 250 µl, 500 µl and 1000 µl
- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Certified RNase, DNase, DNA and Pyrogen free
- Refills dramatically reduce plastic waste compared to pre-racked tips
- Refill top and bottom shells made from 100% recycled plastic
- All components, excluding tips and shrink wrap, are recyclable
- Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

## Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	S1112-1720-C
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	100 - 1,250 µl
Volume	1250 µl
Pack Size	9,600 Tips (100 Refills × 96 Tips)
Color Code Wafer	Blue
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural
Sustainability attributes	Refills dramatically reduce plastic waste compared to pre-racked tips. Refill top and bottom shells made from 100% recycled plastic. All components, excluding tips and shrink wrap wrap, are recyclable. Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



For use with brand + product (SL)

**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima

**Biohit:** ePet, mLine, Proline

**BRAND:** Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

**Capp:** CappAero Comfort, CappAero Universal, Capptronic Single, ecopipette, Singlechannel

**Eppendorf:** Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer

**Gilson:** Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

**Hamilton:** SoftGrip Pipette, SofTouch

**HTL:** Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

**Nichiryo:** 31N, 50F, 50V, NLE, NPLO, NPX

**Rainin:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

**Sartorius:** ePet, mLine, Proline, Proline Plus, Tacta

**Socorex:** Acura 815 manual, Acura 825, Acura electro 925, Calibra digital

**STARLAB:** ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

**Thermo Fisher Scientific Inc.:** Finn Focus, Finnipette Digital, Finnipette F1/F2, Finnipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile), Case

High quality, reliable [pipette tips](#) with a universal tip fit – [TipOne](#)® has it all! If you are looking for a sustainable [pipette](#) tip system with accuracy and precision you can trust, [TipOne](#)® is the tip for you.

[TipOne](#)® refills are a quick and easy way to obtain economical racked tips, and they save on plastic use!

Reduce single-use plastic!













Compared to pre-racked tips, using [TipOne](#)® refills dramatically reduces the amount of single-use plastic going through your laboratory. Quick and easy to load, unwrap and click in place into an empty [TipOne](#)® rack. The refills are supplied in two separate towers of 5x96 tips, so you can rack up 5 or 10 at a time. Tip information is printed on every wafer, providing full traceability when racked. Once racked, and if needed, the refills can be autoclaved.

Join our journey to sustainable pipetting

The clear, polypropylene system components are made from 100% recycled plastic. [TipOne](#)® tip racks are built to last and can be reused again and again with these cost-effective refills. All components of [TipOne](#)® refills, excluding tips and shrink wrap, are recyclable. Sign up for our [TipOne](#)® recycling service!











## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Racks (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Racks × 96 Tips)	S1112-1820
 <p><a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Racks (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 3,840 Tips (40 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	3,840 Tips (40 Racks × 96 Tips)	S1112-1820-C
 <p><a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile)</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Racks × 96 Tips)	S1112-1830
 <p><a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile), Case</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 3,840 Tips (40 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	3,840 Tips (40 Racks × 96 Tips)	S1112-1830-C
 <p><a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Refills × 96 Tips)	S1112-1720
 <p><a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 9,600 Tips (100 Refills × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	9,600 Tips (100 Refills × 96 Tips)	S1112-1720-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Refills (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1112-1730
 <a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Refills (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 3,840 Tips (40 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	3,840 Tips (40 Refills × 96 Tips)	S1112-1730-C
 <a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Bag (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	S1112-1020
 <a href="#">1,250 µl XL graduated TipOne® Tip, Natural, Bags (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 16,000 Tips (16 Bags × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	16,000 Tips (16 Bags × 1000 Tips)	S1112-1020-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)