



0.2 ml PCR Tube, Flat Cap, Natural



Product Highlights

- PCR Tubes manufactured from medical grade polypropylene
- Positive click on closing gives assurance that tubes are fully closed
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Certified free of RNase, Dnase, DNA & Pyrogen-free
- PCR inhibitor free
- Supplied in resealable bags

Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	I1402-8100
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
Pack Size	1,000 Tubes (1 Bag × 1000 Tubes)
Cap style	Flat cap
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Suitable for PCR	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 0.2 ml PCR Tube, Flat Cap, Natural











Starlab 0.2 ml PCR tubes with attached caps provide optimal sealing which minimises fluid loss through evaporation. The caps are easy to open and close yet ensure a tight fit to reduce sample evaporation.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk











All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.2 ml PCR Tube, Domed Cap, Natural Cap style: Domed Color: Natural Pack Size: 1,000 Tubes (1 Bag × 1000 Tubes)  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1402-4300
 0.2 ml PCR Tube, Domed Cap, Mixed Cap style: Domed Color: Mixed Pack Size: 1,000 Tubes (5 Bags × 200 Tubes)  Create PDF Data Sheet	1,000 Tubes (5 Bags × 200 Tubes)	I1402-4308
 0.2 ml PCR Tube, Flat Cap, Natural Cap style: Flat cap Color: Natural Pack Size: 1,000 Tubes (1 Bag × 1000 Tubes)  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1402-8100
 0.2 ml PCR Tube, Flat Cap, Mixed Cap style: Flat cap Color: Mixed Pack Size: 1,000 Tubes (5 Bags × 200 Tubes)  Create PDF Data Sheet	1,000 Tubes (5 Bags × 200 Tubes)	I1402-8108
 0.2 ml PCR Tube, Flat Cap (Xtra-Clear), Natural Cap style: Flat (Xtra-Clear) Color: Natural Pack Size: 1,000 Tubes (1 Bag × 1000 Tubes)  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1402-8200

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.2 ml PCR Tube, Flat Cap, Natural (Case Size) Cap style: Flat cap Color: Natural Pack Size: 10,000 Tubes (10 Bags × 1000 Tubes)  Create PDF Data Sheet	10,000 Tubes (10 Bags × 1000 Tubes)	I1402-8100-C
 0.2 ml PCR Tube, Flat Cap, Mixed (Case Size) Cap style: Flat cap Color: Mixed Pack Size: 10,000 Tubes (50 Bags × 200 Tubes)  Create PDF Data Sheet	10,000 Tubes (50 Bags × 200 Tubes)	I1402-8108-C
 0.2ml PCR Tube, Domed Shielded Cap, Natural Cap style: Domed shielded cap Color: Natural Pack Size: 1,000 Tubes (1 Bag × 1000 Tubes)  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	S1602-4300
 0.2 ml PCR Tube, Flat Shielded Cap, Natural Cap style: Flat shielded cap Color: Natural Pack Size: 1,000 Tubes (1 Bag × 1000 Tubes)  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	B1402-5500

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk