



## 0.2 ml PCR Tube, Flat Shielded Cap, Natural



### Product Highlights

- PCR Tubes manufactured from medical grade polypropylene
- Positive click on closing gives assurance that tubes are fully closed
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Certified free of RNase, Dnase, DNA & Pyrogen-free
- PCR inhibitor free
- The integral shield on the cap prevents accidental finger contamination with the surface of the lid
- Supplied in resealable bags

## Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                    |                                  |
|--------------------|----------------------------------|
| Art. No.           | B1402-5500                       |
| Pack Size          | 1,000 Tubes (1 Bag x 1000 Tubes) |
| Sterile            | No                               |
| Autoclavable       | Yes                              |
| Volume             | 0.2 ml                           |
| Cap style          | Flat shielded cap                |
| Material           | Polypropylene                    |
| DNA free           | Yes                              |
| DNase free         | Yes                              |
| Free of endotoxins | Yes                              |
| Human DNA free     | Yes                              |
| PCR inhibitor free | Yes                              |
| Pyrogen free       | Yes                              |
| RNase free         | Yes                              |
| Color              | Natural                          |
| Suitable for PCR   | yes                              |


STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 0.2 ml PCR Tube, Flat Shielded Cap, Natural

Starlab 0.2 ml PCR tubes with attached caps provide optimal sealing which minimises fluid loss through evaporation. The caps are easy to open and close yet ensure a tight fit to reduce sample evaporation.




## All Variations

| PRODUCT NAME                                                                                                                                                                                                                                                                                                                                                                                                 | PACKAGING SIZE                   | ART. NO.   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|------------|
| <br><a href="#">0.2 ml PCR Tube, Domed Cap, Natural</a><br><b>Cap style:</b> Domed<br><b>Color:</b> Natural<br><b>Pack Size:</b> 1,000 Tubes<br>(1 Bag × 1000 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                             | 1,000 Tubes (1 Bag × 1000 Tubes) | I1402-4300 |
| <br><a href="#">0.2 ml PCR Tube, Domed Cap, Mixed</a><br><b>Cap style:</b> Domed<br><b>Color:</b> Mixed<br><b>Pack Size:</b> 1,000 Tubes<br>(5 Bags × 200 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                                 | 1,000 Tubes (5 Bags × 200 Tubes) | I1402-4308 |
| <br><a href="#">0.2 ml PCR Tube, Flat Cap, Natural</a><br><b>Cap style:</b> Flat cap<br><b>Color:</b> Natural<br><b>Pack Size:</b> 1,000 Tubes<br>(1 Bag × 1000 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                        | 1,000 Tubes (1 Bag × 1000 Tubes) | I1402-8100 |
| <br><a href="#">0.2 ml PCR Tube, Flat Cap, Mixed</a><br><b>Cap style:</b> Flat cap<br><b>Color:</b> Mixed<br><b>Pack Size:</b> 1,000 Tubes<br>(5 Bags × 200 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                           | 1,000 Tubes (5 Bags × 200 Tubes) | I1402-8108 |
| <br><a href="#">0.2 ml PCR Tube, Flat Cap (Xtra-Clear), Natural</a><br><b>Cap style:</b> Flat (Xtra-Clear)<br><b>Color:</b> Natural<br><b>Pack Size:</b> 1,000 Tubes<br>(1 Bag × 1000 Tubes)<br> <a href="#">Create PDF Data Sheet</a> | 1,000 Tubes (1 Bag × 1000 Tubes) | I1402-8200 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



| PRODUCT NAME                                                                        |                                                                                                                                                                                                                                                                                                                        | PACKAGING SIZE                      | ART. NO.     |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|--------------|
|    | <a href="#">0.2 ml PCR Tube, Flat Cap, Natural (Case Size)</a><br><b>Cap style:</b> Flat cap<br><b>Color:</b> Natural<br><b>Pack Size:</b> 10,000 Tubes<br>(10 Bags × 1000 Tubes)<br><br> <a href="#">Create PDF Data Sheet</a>       | 10,000 Tubes (10 Bags × 1000 Tubes) | I1402-8100-C |
|    | <a href="#">0.2 ml PCR Tube, Flat Cap, Mixed (Case Size)</a><br><b>Cap style:</b> Flat cap<br><b>Color:</b> Mixed<br><b>Pack Size:</b> 10,000 Tubes<br>(50 Bags × 200 Tubes)<br><br> <a href="#">Create PDF Data Sheet</a>            | 10,000 Tubes (50 Bags × 200 Tubes)  | I1402-8108-C |
|   | <a href="#">0.2ml PCR Tube, Domed Shielded Cap, Natural</a><br><b>Cap style:</b> Domed shielded cap<br><b>Color:</b> Natural<br><b>Pack Size:</b> 1,000 Tubes<br>(1 Bag × 1000 Tubes)<br><br> <a href="#">Create PDF Data Sheet</a> | 1,000 Tubes (1 Bag × 1000 Tubes)    | S1602-4300   |
|  | <a href="#">0.2 ml PCR Tube, Flat Shielded Cap, Natural</a><br><b>Cap style:</b> Flat shielded cap<br><b>Color:</b> Natural<br><b>Pack Size:</b> 1,000 Tubes<br>(1 Bag × 1000 Tubes)<br><br> <a href="#">Create PDF Data Sheet</a>  | 1,000 Tubes (1 Bag × 1000 Tubes)    | B1402-5500   |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)