



0.2ml 8-Strip PCR Tubes, Individually Attached Domed Shielded Caps



Product Highlights

- Manufactured from medical grade polypropylene
- Attached caps eliminates the possibility of cross-contamination through accidental tube–cap mix up
- PCR inhibitor free
- Certified RNase, Dnase, DNA and Pyrogen-free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	S1602-2900
Pack Size	120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
No. of tubes	8
Cap style	Individually attached domed shielded caps
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Suitable for PCR	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 0.2ml 8-Strip PCR Tubes, Individually Attached Domed Shielded Caps









Starlab PCR Strips are manufactured from medical grade polypropylene. The thin-walled tubes ensure efficient and homogenous heat transfer. The attached 'crocodile' style strip caps eliminate the possibility of cross-contamination from cap mix ups. Available in 8- or 12-Strip Tubes and Caps.

Low-Profile Strips (0.1 ml): 15.24 mm high, for use in cyclers with low-profile capability only.

Standard Strips (0.2 ml): 20.83 mm high, for use in the majority of cyclers.













All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear) Cap style: Individually attached flat caps (X-Clear) Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)	A1402-3700
 0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear)(Case Size) Cap style: Individually attached flat caps (X-Clear) Volume: 0.2 ml Pack Size: 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)  Create PDF Data Sheet	1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)	A1402-3700-C
 0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural Cap style: Individually attached domed caps Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)	I1402-2900
 8-Strip PCR Tubes with Individual Attached Caps, 0.2 ml (mixed) Cap style: Individually attached domed caps Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)	I1402-2908

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.1 ml 8-Strip PCR Tubes, Low-Profile, Individually Attached Flat Caps (Xtra-Clear) Cap style: Individually attached flat caps – (Xtra-Clear), low-profile Volume: 0.1 ml Pack Size: 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)	A1402-3800
 0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural (Case Size) Cap style: Individually attached domed caps Volume: 0.2 ml Pack Size: 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)  Create PDF Data Sheet	1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)	I1402-2900-C
 0.2ml 8-Strip PCR Tubes, Individually Attached Domed Shielded Caps Cap style: Individually attached domed shielded caps Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)	S1602-2900
 0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip Cap style: Attached flat cap strip (Xtra-Clear) Volume: 0.2 ml Pack Size: 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)  Create PDF Data Sheet	125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)	I1402-2700
 0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip Cap style: Attached domed cap strip Volume: 0.2 ml Pack Size: 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)  Create PDF Data Sheet	125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)	I1402-2800

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME		PACKAGING SIZE	ART. NO.
	0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip (Case Size)	1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)	I1402-2700-C
	Cap style: Attached flat cap strip (Xtra-Clear) Volume: 0.2 ml Pack Size: 1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)  Create PDF Data Sheet		
	0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip (Case Size)	1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)	I1402-2800-C
	Cap style: Attached domed cap strip Volume: 0.2 ml Pack Size: 1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)  Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk