



## 0.2 ml 8-Strip PCR Tubes, White with Attached Flat (Xtra-Clear) Cap Strip



### Product Highlights

- Manufactured from medical grade polypropylene
- Individually attached caps eliminates cross-contamination from cap mix-up
- PCR inhibitor free
- Suitable for Realtime PCR (quantitative PCR = qPCR)
- Certified RNase, Dnase, DNA and Pyrogen-free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

## Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	I1409-2700
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
Pack Size	125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)
No. of tubes	8
Cap style	Attached flat cap strip (Xtra-Clear)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	White
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 0.2 ml 8-Strip PCR Tubes, White with Attached Flat (Xtra-Clear) Cap Strip

### Help eliminate cross-contamination with PCR Strips with Attached Cap Strip

Discover convenience and efficiency with Starlab's PCR strips with attached cap strip, designed for mid-throughput settings. These PCR strips ensure all tubes are equally protected, supporting consistent results across samples, as well as simplifying handling to increase workflow efficiency. The attached 'crocodile' style strip caps eliminate the possibility of cross-contamination from cap mix ups.

The opaque white tubes maximise reflectance for enhanced signal, and the ultra-thin, uniform walls ensure optimal heat transfer. Combined with the Xtra-Clear cap strip for optimal signal transmission and detection, these strips are ideal for qPCR.

0.2ml tubes are suitable for use in the majority of cyclers. Certified free of RNase, DNase, DNA and PCR inhibitors, and tested Pyrogen free.