



96-Well PCR Plate, Semi-Skirted, for Roche® Lightcycler®



Product Highlights

- Low-profile, semi-skirted PCR plate, with a maximum capacity of 200 μ l
- Opaque white plate ideal for qPCR. White plates provide an optimal signal-to-noise ratio for fluorescence-based assays.
- For use with Roche® Lightcycler®

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|----------------------|--------------------------------------|
| Art. No. | I1402-9909 |
| Sterile | No |
| Autoclavable | No |
| Volume | 200 µl |
| Pack Size | 10 Plates (1 Box × 10 Plates) |
| Number of wells | 96 wells |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| PCR inhibitor free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Color | White |
| Cut corner | H12 |
| Elevated wells | no |
| Frame Design | Semi-skirted |
| Matrix | Printed |
| Profile | Low profile, for Roche® Lightcycler® |
| Raised rim | no |
| Suitable for PCR | yes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



| | |
|-----------------------------------|-----|
| Suitable for Real-Time PCR (qPCR) | yes |
|-----------------------------------|-----|

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 96-Well PCR Plate, Semi-Skirted, for Roche® Lightcycler®

Starlab's semi-skirted PCR plates are made from medical polypropylene to provide the quality assurance you need. More rigid than non-skirted plates, they also offer better handling convenience. The plates have clearly marked wells to provide easy identification and sample tracking, as well as compatibility with a wide range of thermal cyclers.

Supplied in convenient, resealable bags to minimise contamination risk. Certified free of RNase, DNase, DNA and PCR Inhibitors. Tested Pyrogen free.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



Accessories

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|-----------------------------|------------|
|  Aluminium StarSeal (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-9792 |
|  Aluminium Sealing Film, 60 µm (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-0792 |
|  Clear Polyolefin StarSeal (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-9793 |
|  Clear Polypropylene Seal (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-0793 |
|  Xtra-Clear Advanced Polyolefin StarSeal (qPCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-9795 |
|  Polyester PCR Sealing Film, Clear | 100 Pcs. (1 Box × 100 Pcs.) | E2796-0100 |
|  Polypropylene PCR Sealing Film Strips, Clear | 200 Pcs. (1 Box × 200 Pcs.) | E2796-2850 |
|  Self-adhesive sealing films & foils for PCR | 100 Pcs. (1 Box × 100 Pcs.) | E2796-1270 |
|  Silicone Sealing Mat for 96-Well PCR Plate | 5 Pcs. (1 Box × 5 Pcs.) | E1403-0000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk