



96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white (Case Size)



Product Highlights

- Approximate maximum capacity for all 96-well plates is 350 µl for standard height plates, or 200 µl for low-profile plates
- Available as natural, white or black plates. White plates are ideal for qPCR, giving optimal signal-to-noise ratio for fluorescent based assays. Black plates help to minimize light diffusion.
- Specialist plates (for FAST® Systems, qPCR)

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	E1403-7709-C
Sterile	No
Autoclavable	Yes
Volume	200 µl
Pack Size	100 Plates (10 Boxes × 10 Plates)
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	White
Cut corner	A1
Elevated wells	no
Frame Design	Semi-skirted
Matrix	Printed
Profile	Low (for FAST® Systems), raised rim
Raised rim	yes
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white (Case Size)

STARLAB Semi-Skirted PCR Plates feature a 7 mm high skirt around the complete plate edge.

The skirt improves the stability of the plate and makes them more suited e. g. for automated applications. Furthermore a labelling on the plate edge e. g. by barcode-labels is possible. Plates with elevated wells reduce the contamination risk by cross-contamination.









- > Autoclavable
- > Supplied in resealable bags
- > Convenient small pack sizes
- > Certified RNase, DNase, DNA and PCR Inhibitor free
- > Endotoxin free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk







All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), natural</p> <p>Profile: Low (for FAST® Systems), raised rim Color: Natural Pack Size: 10 Plates (1 Box × 10 Plates)</p> <p> Create PDF Data Sheet</p>	10 Plates (1 Box × 10 Plates)	E1403-7700
 <p>96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), natural (Case Size)</p> <p>Profile: Low (for FAST® Systems), raised rim Color: Natural Pack Size: 100 Plates (10 Boxes × 10 Plates)</p> <p> Create PDF Data Sheet</p>	100 Plates (10 Boxes × 10 Plates)	E1403-7700-C
 <p>96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white</p> <p>Profile: Low (for FAST® Systems), raised rim Color: White Pack Size: 10 Plates (1 Box × 10 Plates)</p> <p> Create PDF Data Sheet</p>	10 Plates (1 Box × 10 Plates)	E1403-7709
 <p>96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white (Case Size)</p> <p>Profile: Low (for FAST® Systems), raised rim Color: White Pack Size: 100 Plates (10 Boxes × 10 Plates)</p> <p> Create PDF Data Sheet</p>	100 Plates (10 Boxes × 10 Plates)	E1403-7709-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk











PRODUCT NAME	PACKAGING SIZE	ART. NO.
 96-Well PCR Plate, Semi-Skirted, Raised Rim, natural Profile: Raised rim Color: Natural Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E1403-8200
 96-Well PCR Plate, Semi-Skirted, Raised Rim, natural (Case Size) Profile: Raised rim Color: Natural Pack Size: 100 Plates (10 Boxes × 10 Plates)  Create PDF Data Sheet	100 Plates (10 Boxes × 10 Plates)	E1403-8200-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Silicone Sealing Mat for 96-Well PCR Plate	5 Pcs. (1 Box × 5 Pcs.)	E1403-0000
 Aluminium StarSeal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9792
 Aluminium Sealing Film, 60 µm (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-0792
 Clear Polyolefin StarSeal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9793
 Clear Polypropylene Seal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-0793
 Xtra-Clear Advanced Polyolefin StarSeal (qPCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9795
 Polyester PCR Sealing Film, Clear	100 Pcs. (1 Box × 100 Pcs.)	E2796-0100
 Polypropylene PCR Sealing Film Strips, Clear	200 Pcs. (1 Box × 200 Pcs.)	E2796-2850
 Self-adhesive sealing films & foils for PCR	100 Pcs. (1 Box × 100 Pcs.)	E2796-1270

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk