



96-Well PCR Plate, Skirted, Low Profile, White



Product Highlights

- Approximate maximum capacity for all 96-well plates is 350 μ l for standard height plates, or 200 μ l for low-profile plates
- Holes on sides for robotics handling
- Highly-visible printed matrix
- Certified RNase, DNase, DNA and Pyrogen-free
- Opaque white plates available for qPCR
- PCR Inhibitor free
- –
- –

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E1403-5209
Pack Size	10 Plates (1 Pack × 10 Plates)
Sterile	No
Autoclavable	Yes
Volume	200 µl
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	White
Cut corner	A12
Frame Design	Skirted
Matrix	Printed
Profile	Low
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk












More informations about 96-Well PCR Plate, Skirted, Low Profile, White

Starlab Skirted 96- and 384-Well PCR Plates feature a 15 mm high full skirt around the complete plate edge. The skirt provides an extremely high plate stability, which is beneficial specifically in automated "High-Throughput" applications. These type of plates offer the largest area for marking on edge of the plates e. g. by barcodes. Starlab Skirted PCR Plates are available in standard or low profile format.



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Silicone Sealing Mat for 96-Well PCR Plate	5 Pcs. (1 Box × 5 Pcs.)	E1403-0000
 Aluminium StarSeal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9792
 Aluminium Sealing Film, 60 µm (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-0792
 Clear Polyolefin StarSeal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9793
 Clear Polypropylene Seal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-0793
 Xtra-Clear Advanced Polyolefin StarSeal (qPCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9795
 Polyester PCR Sealing Film, Clear	100 Pcs. (1 Box × 100 Pcs.)	E2796-0100
 Polypropylene PCR Sealing Film Strips, Clear	200 Pcs. (1 Box × 200 Pcs.)	E2796-2850
 Self-adhesive sealing films & foils for PCR	100 Pcs. (1 Box × 100 Pcs.)	E2796-1270

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk