



## 96-Well PCR Plate, Skirted, Low Profile (Case Size)



### Product Highlights

- Approximate maximum capacity for all 96-well plates is 350  $\mu$ l for standard height plates, or 200  $\mu$ l for low-profile plates
- Holes on sides for robotics handling
- Highly-visible printed matrix
- Certified RNase, DNase, DNA and Pyrogen-free
- PCR Inhibitor free
- -
- -

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	E1403-5200-C
Sterile	No
Autoclavable	Yes
Volume	200 µl
Pack Size	100 Plates (10 Packs × 10 Plates)
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Cut corner	A12
Frame Design	Skirted
Matrix	Printed
Profile	Low
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 96-Well PCR Plate, Skirted, Low Profile (Case Size)

Starlab Skirted 96- and 384-Well [PCR Plates](#) feature a 15 mm high full skirt around the complete plate edge. The skirt provides an extremely high plate stability, which is beneficial specifically in automated "High-Throughput" applications. These type of plates offer the largest area for marking on edge of the plates e. g. by barcodes. Starlab Skirted [PCR Plates](#) are available in standard or low profile format.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Aluminium StarSeal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9792
 <a href="#">Aluminium Sealing Film, 60 µm (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0792
 <a href="#">Clear Polyolefin StarSeal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9793
 <a href="#">Clear Polypropylene Seal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0793
 <a href="#">Xtra-Clear Advanced Polyolefin StarSeal (qPCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9795
 <a href="#">Polyester PCR Sealing Film, Clear</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0100
 <a href="#">Polypropylene PCR Sealing Film Strips, Clear</a>	200 Pcs. (1 Box × 200 Pcs.)	E2796-2850
 <a href="#">Self-adhesive sealing films &amp; foils for PCR</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-1270
 <a href="#">Silicone Sealing Mat for 96-Well PCR Plate</a>	5 Pcs. (1 Box × 5 Pcs.)	E1403-0000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)