



## 96-Well PCR Plate, Semi-Skirted, Straight Edges, with Barcode



### Product Highlights

- Code 128 Barcoding plus human-readable number
- All barcodes seven digits long
- Low, minimum ordering quantity of 100 plates
- Certified RNase, DNase, DNA and Pyrogen-free
- PCR Inhibitor free

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	I1402-9800-BC
Sterile	No
Autoclavable	Yes
Volume	350 µl
Pack Size	100 Plates (10 Packs × 10 Plates)
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Cut corner	A12
Elevated wells	no
Frame Design	Semi-skirted
Profile	Standard height, straight edges
Raised rim	no
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)













## More informations about 96-Well PCR Plate, Semi-Skirted, Straight Edges, with Barcode

### Customised barcoding is also available

If you require a specific sequence start to your barcoding or require barcoding on plates other than those shown above, we can provide you with a quote for customised barcoding. Note that the minimum order quantity for custom barcoding is 2000 plates. Contact Starlab or your Starlab Account Manager for more details.



## All Variations










PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">384-Well PCR Plate for Roche Lightcycler, Skirted, with Barcode</a> <b>Profile:</b> For Roche® Lightcycler <b>Frame Design:</b> Skirted  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Packs × 10 Plates)	E1042-9900-BC
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, with Barcode</a> <b>Profile:</b> Raised rim <b>Frame Design:</b> Semi-skirted  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Packs × 10 Plates)	E1403-8200-BC
 <a href="#">96-Well PCR Plate, White, Semi-Skirted, for Roche® Lightcycler, with Barcode</a> <b>Profile:</b> Low profile, for Roche® Lightcycler® <b>Frame Design:</b> Semi-skirted  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Packs × 10 Plates)	I1402-9909-BC
 <a href="#">96-Well PCR Plate, Skirted, Low-Profile, with Barcode</a> <b>Profile:</b> Low <b>Frame Design:</b> Skirted  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Packs × 10 Plates)	E1403-5200-BC
 <a href="#">96-Well PCR Plate, Semi-Skirted, Straight Edges, with Barcode</a> <b>Profile:</b> Standard height, straight edges <b>Frame Design:</b> Semi-skirted  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Packs × 10 Plates)	I1402-9800-BC

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Aluminium Sealing Film, 60 µm (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0792
 <a href="#">Clear Polyolefin StarSeal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9793
 <a href="#">Clear Polypropylene Seal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0793
 <a href="#">Xtra-Clear Advanced Polyolefin StarSeal (qPCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9795
 <a href="#">Polyester PCR Sealing Film, Clear</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0100
 <a href="#">Polypropylene PCR Sealing Film Strips, Clear</a>	200 Pcs. (1 Box × 200 Pcs.)	E2796-2850
 <a href="#">Self-adhesive sealing films &amp; foils for PCR</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-1270
 <a href="#">Aluminium StarSeal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9792
 <a href="#">Silicone Sealing Mat for 96-Well PCR Plate</a>	5 Pcs. (1 Box × 5 Pcs.)	E1403-0000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)