




12-Strip Caps, Domed



Product Highlights

- Manufactured from medical grade polypropylene
- 12-Strip Caps can be easily split into smaller sections by snapping apart
- Certified RNase, DNase, DNA & Pyrogen-free
- PCR inhibitor free
- Supplied in resealable bags

Product Applications

PCR

Realtime PCR (quantitative PCR = qPCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I1400-1200
Pack Size	80 Pcs. (1 Bag × 80 Pcs.)
Sterile	No
Autoclavable	Yes
Cap style	Domed
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Suitable for PCR	yes
For use with	PCR Tubes & Plates

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 12-Strip Caps, Domed









12-strip PCR caps. A range of different types of PCR strip tubes are available separately, so you can mix and match tube and cap style to suit your needs.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 12-Strip Caps, Domed Cap style: Domed Color: Natural Pack Size: 80 Pcs. (1 Bag × 80 Pcs.)  Create PDF Data Sheet	80 Pcs. (1 Bag × 80 Pcs.)	I1400-1200
 12-Strip PCR Caps, Domed, Natural (Case Size) Cap style: Domed Color: Natural Pack Size: 800 Pcs. (10 Bags × 80 Pcs.)  Create PDF Data Sheet	800 Pcs. (10 Bags × 80 Pcs.)	I1400-1200-C
 12-Strip PCR Caps, Domed, Mixed Cap style: Domed Color: Mixed Pack Size: 80 Pcs. (5 Bags × 16 Pcs.)  Create PDF Data Sheet	80 Pcs. (5 Bags × 16 Pcs.)	I1400-1208
 12-Strip PCR Caps, Flat (Xtra-Clear), Natural Cap style: Flat (Xtra-Clear) Color: Natural Pack Size: 80 Pcs. (1 Bag × 80 Pcs.)  Create PDF Data Sheet	80 Pcs. (1 Bag × 80 Pcs.)	I1400-1290

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk