



2.0 ml 96 Deepwell Plate, Round Wells with Round Bottoms (Sterile)



Product Highlights

- Virgin polypropylene plates with SBS footprint
- Alpha-numeric reference grid (moulded)
- Plate temperature use: -80 °C to 121 °C.
- Autoclavable
- Max. centrifugation rate: 5000 x g

Product Applications

Sample preparation and storage
DNA and RNA isolation with high throughput
Plasmid purification with high throughput
Storage of cDNA or genomic banks
Storage of siRNA or oligonucleotide libraries
Protein analysis
Bacteria and yeast cultivation
Dilutions series

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2896-2110
Pack Size	10 Plates (2 Packs × 5 Plates)
Sterile	Yes
Autoclavable	Yes
Volume	2.0 ml
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
RNase free	Yes
Max. temp.	121 °C
Temperature range	121 °C
Cut corner	H1/H12
g-Force max.	5,000 × g
Profile	Round wells with round bottoms

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 2.0 ml 96 Deepwell Plate, Round Wells with Round Bottoms (Sterile)


Virgin polypropylene plates with SBS footprint. Alpha-numeric reference grid (moulded).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 2.5ml 96 Deepwell Plate, Square Wells with Round Bottoms Volume: 2.5 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet	20 Plates (1 Box × 20 Plates)	E2896-2500
 10 ml 24 Deepwell Plate, Square Wells with V-Shape Bottoms Volume: 10 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet	20 Plates (1 Box × 20 Plates)	E2824-0000
 5.0 ml 48 Deepwell Plate, Rectangle Wells with V-Shape Bottoms Volume: 5.0 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet	20 Plates (1 Box × 20 Plates)	E2848-0000
 1.2 ml 96 Deepwell Plate, Round Wells with Round Bottoms, Low-Profile Volume: 1.2 ml Sterile: No Pack Size: 50 Plates (1 Box × 50 Plates)  Create PDF Data Sheet	50 Plates (1 Box × 50 Plates)	E2896-0120
 2.2 ml 96 Deepwell Plate, Square Wells with V-Shaped Bottoms Deepwell Plates Volume: 2.20 ml Sterile: No Pack Size: 50 Plates (1 Box × 50 Plates)  Create PDF Data Sheet	50 Plates (1 Box × 50 Plates)	E2896-0220

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




PRODUCT NAME	PACKAGING SIZE	ART. NO.
 2.0 ml 96 Deepwell Plate, Round Wells with Round Bottoms (Sterile) Volume: 2.0 ml Sterile: Yes Pack Size: 10 Plates (2 Packs × 5 Plates)  Create PDF Data Sheet	10 Plates (2 Packs × 5 Plates)	E2896-2110
 2.0 ml 96 Deepwell Plate, Round Wells with Round Bottoms Volume: 2.0 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2896-2100
 1.2 ml 96 Deepwell Plate, Square Wells with Conical Bottoms Volume: 1.2 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2896-1820
 2.2 ml 96 Deepwell Plate, Square Wells with Conical Bottoms Volume: 2.2 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2896-1810
 0.6 ml 96 Deepwell Plate, Round Wells with Conical Bottoms Volume: 0.6 ml Sterile: No Pack Size: 5 Plates (1 Box × 5 Plates)  Create PDF Data Sheet	5 Plates (1 Box × 5 Plates)	E2896-0600
 1.1 ml 96 Deepwell Plate, Round Wells with Round Bottoms, Low-Profile (Sterile) Volume: 1.1 ml Sterile: Yes Pack Size: 10 Plates (2 Packs × 5 Plates)  Create PDF Data Sheet	10 Plates (2 Packs × 5 Plates)	E2896-0110

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.


Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 240 µl 384 Deepwell Plate, Square Wells with V-Shape Bottoms Volume: 0.24 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2384-2400

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Silicone Sealing Mat for 96-Well Deepwell-Plates, 1 ml / 2 ml, with Round Wells, piercable, self-sealing	10 Pcs. (1 Box × 10 Pcs.)	E2896-1100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk