



200 µl UltraPoint® graduated TipOne® Tip, Natural, Stack racks (non-sterile), Case



Product Highlights

- **Graduated** at 10 µl, 50 µl and 100 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified** RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	S1113-1200-C
Pack Size	9,600 Tips (10 Stack racks × 960 Tips)
Color Code Wafer	Yellow
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Tip Type	UltraPoint®
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with	<p>Abimed: BioPette-E, Discovery Comfort, LABMATE Optima</p> <p>Biohit: eLine, ePet, mLine, Picus, Proline</p> <p>BRAND: Transferpette, Transferpette Electronic, Transferpette S</p> <p>Capp: CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel</p> <p>Eppendorf: Reference, Reference 2, Research, Research plus, Research pro, Xplorer</p> <p>Gilson: Pipetman Concept, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p>Hamilton: SoftGrip Pipette, SofTouch</p> <p>HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p>Nichiryo: 00-NP7, 50F, 50V, NLE, NPLO, NPX</p> <p>Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p>Sartorius: eLine, ePet, mLine, Picus, Picus Nxt, Proline, Proline Plus, Tacta</p> <p>Socorex: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital</p> <p>STARLAB: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p>Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finn Focus Univ., Finn timer Digital, Finn timer F1/F2, Finn timer Novus, Finn timer Novus Turquoise, Finn timer Novus Yellow</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 200 µl UltraPoint® graduated TipOne® Tip, Natural, Stack racks (non-sterile), Case













UltraPoint orifice

Thin-wall orifice helps prevent sample hang up and makes it easier to touch off the last drop of sample

The compact Stack Rack provides 960 tips in one convenient unit, making an excellent choice for high throughput or multi-channel use. The wide base provides stability during use and the transparent cover protects the tips while not in use. The Stack Rack is autoclavable and all components are made from polypropylene and can be recycled.















All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Rack (non-sterile) Sterile: No Pack Size: 960 Tips (10 Rack × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Rack × 96 Tips)	S1113-1800
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1113-1800-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1113-1810
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1113-1810-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1113-1700
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1113-1700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl UltraPoint® graduated TipOne® Tip, Natural Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1113-1710
 200 µl UltraPoint® graduated TipOne® Tip, Natural Refills (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Refills × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Refills × 96 Tips)	S1113-1710-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)  Create PDF Data Sheet	960 Tips (1 Stack rack × 960 Tips)	S1113-1200
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)  Create PDF Data Sheet	9,600 Tips (10 Stack racks × 960 Tips)	S1113-1200-C
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	S1113-1000
 200 µl UltraPoint® graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 25,000 Tips (25 Bags × 1000 Tips)  Create PDF Data Sheet	25,000 Tips (25 Bags × 1000 Tips)	S1113-1000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk