



StarTub PS Reagent Reservoir, 25 ml, Individually Wrapped (Sterile)



Product Highlights

- White polystyrene (PS) reservoir
- Sterile
- For Single use
- Capacity: 25 ml
- free of detectable DNA, DNase and RNase and Endotoxins

Product Applications

Reagent takeup with multichannel pipette

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2310-1025
Dimensions (W × D × H)	145 mm × 70 mm × 26 mm
Dimensions (W × D)	145 mm × 70 mm
Sterile	Yes
Autoclavable	No
Volume	25 ml
Pack Size	50 Pcs. (1 Box × 50 Pcs.)
Material	Polystyrene
Color	White
Bottoms	Conical
Reusable	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about StarTub PS Reagent Reservoir, 25 ml, Individually Wrapped (Sterile)









25 ml polystyrene, sterile individually wrapped reservoirs with V-shaped bottom. Graduations.
Dimensions (WxDxH): 145 x 70 x 26 mm.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 StarTub PS Reagent Reservoir, 55 ml, Individually Wrapped (Sterile) Sterile: Yes Volume: 55 ml  Create PDF Data Sheet	50 Pcs. (1 Box × 50 Pcs.)	E2310-1010
 StarTub PS Reagent Reservoir, 100 ml, Bulk Sterile: No Volume: 100 ml  Create PDF Data Sheet	100 Pcs. (1 Box × 100 Pcs.)	E2310-2000
 StarTub PS Reagent Reservoir, 25 ml, Individually Wrapped (Sterile) Sterile: Yes Volume: 25 ml  Create PDF Data Sheet	50 Pcs. (1 Box × 50 Pcs.)	E2310-1025
 StarTub PS Reagent Reservoir, 100 ml, Individually Wrapped (Sterile) Sterile: Yes Volume: 100 ml  Create PDF Data Sheet	50 Pcs. (1 Box × 50 Pcs.)	E2310-1100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk