



CAPP® Origami Reagent Reservoir, 40ml



Product Highlights

- Flat-packed durable polypropylene reservoirs
- V-shaped bottom helps prevent bubbles forming

Product Applications

Reagent takeup with multichannel pipette

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2310-0510
Pack Size	50 Pcs. (1 Box × 50 Pcs.)
Volume	40 ml
Material	Polypropylene
Bottoms	Conical

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about **CAPP® Origami Reagent Reservoir, 40ml**

Flat-packed durable polypropylene reservoirs ready for you to flick into shape.





The V-shaped bottom helps prevent bubbles forming and save on reagents.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	CAPP® Origami Reagent Reservoir, 40ml Volume: 40 ml	50 Pcs. (1 Box × 50 Pcs.)	E2310-0510
	Create PDF Data Sheet		
	CAPP® Origami Reagent Reservoir, 30ml Volume: 30 ml	50 Pcs. (1 Box × 50 Pcs.)	E2310-0505
	Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk