




Lid for 96-Well PS Microplates, non-sterile



Product Highlights

- Crystal-clear polystyrene lids for 96-well PS microplates
- For use with plates: E2996-1600, E2996-1700 and E2996-1800
- Aseptically produced in a fully automated and bio-controlled environment
- Supplied in convenient mini-packs of 20 lids

Product Applications

Sample preparation
Homogenous standard assays
Extensive screening assays
ELISA

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2996-1900
Pack Size	100 Lids (5 Packs × 20 Lids)
Sterile	No
Autoclavable	No
Material	Polystyrene
For use with	Plates: E2996-1600, E2996-1700 and E2996-1800

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk







More informations about Lid for 96-Well PS Microplates, non-sterile

Aseptically* produced lids made from crystal clear polystyrene (PS). For use with plates E2996-1600, E2996-1700 and E2996-1800 only. Supplied in convenient mini-packs of 20 lids.

*Aseptically produced in a fully automated and bio-controlled environment.



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	PS Lid for 96-Well PS Microplates, sterile, Ind. Wrapped Sterile: Yes	50 Individually wrapped sachets (1 Pack × 50 Individually wrapped sachets)	E2996-1910
	 Create PDF Data Sheet		
	Lid for 96-Well PS Microplates, non-sterile Sterile: No	100 Lids (5 Packs × 20 Lids)	E2996-1900
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk