



## 200 µl Bevelled TipOne® Tip, Natural, Stack rack (sterile)



### Product Highlights

- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and certified RNase, DNase, DNA and Pyrogen free
- All components, excluding tips and shrink wrap, are recyclable

## Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	S1111-1210
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Pack Size	960 Tips (1 Stack rack × 960 Tips)
Color Code Wafer	Yellow
Tip Type	bevelled
Material	Polypropylene
Graduations	No
Color	Natural
Sustainability attributes	All components, excluding tips and shrink wrap, are recyclable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



For use with brand + product (SL)

**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima

**Biohit:** eLine, ePet, mLine, Picus, Proline

**BRAND:** Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

**Capp:** CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel

**Eppendorf:** Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer

**Gilson:** Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

**Hamilton:** SoftGrip Pipette, SofTouch

**HTL:** Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

**Nichiryo:** 00-NP7, 50F, 50V, NLE, NPLO, NPX

**Rainin:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

**Sartorius:** eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

**Socorex:** Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

**STARLAB:** ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

**Thermo Fisher Scientific Inc.:** Finn Focus, Finn Focus Micro, Finn Focus Univ., Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Turquoise, Finnpipette Novus Yellow

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 200 µl Bevelled TipOne® Tip, Natural, Stack rack (sterile)

High quality, reliable [pipette tips](#) with a universal tip fit – [TipOne®](#) has it all! If you are looking for a sustainable [pipette](#) tip system with accuracy and precision you can trust, [TipOne®](#) is the one tip for you.

Sterile [TipOne®](#) tips are validated according to EN ISO 11137-2:2007.

960 tips in one, convenient and compact unit!













The Stack Rack is ideal for high-throughput, multi-channel use or space-limited lab environments. Each layer is sealed and protected until it is needed. The Stack Rack is designed to provide quick and easy access to the next layer of tips, without the need to stop pipetting to reload another rack. The Stack Rack drastically reduces plastic use compared to racked tips.

Join our journey to sustainable pipetting

When you have used all your tips, the coloured wafers, base, layers and lid can be recycled. Sign up for our [TipOne®](#) recycling service!















## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Racks × 96 Tips)	S1111-1800
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	7,680 Tips (80 Racks × 96 Tips)	S1111-1800-C
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (sterile)</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Racks × 96 Tips)	S1111-1810
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (sterile), Case</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	7,680 Tips (80 Racks × 96 Tips)	S1111-1810-C
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Refills × 96 Tips)	S1111-1700
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 9,600 Tips (100 Refills × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	9,600 Tips (100 Refills × 96 Tips)	S1111-1700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack rack (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (1 Stack rack × 960 Tips)	S1111-1200
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack racks (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-1200-C
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack rack (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (1 Stack rack × 960 Tips)	S1111-1210
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack racks (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-1210-C
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Bag (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	S1111-1000
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Bags (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 25,000 Tips (25 Bags × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	25,000 Tips (25 Bags × 1000 Tips)	S1111-1000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl Bevelled TipOne® Tip, Natural Refills (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1111-1710
 <a href="#">200 µl Bevelled TipOne® Tip, Natural Refills (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Refills × 96 Tips)	S1111-1710-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)