




## 1.0 ml Graduated Pasteur Pipette



### Product Highlights

- Low affinity surface
- Attached bulb avoids cross-contamination
- One-step biohazard disposal
- Shatter-proof low-density polyethylene, which is non-toxic and inert
- Can be frozen in liquid nitrogen
- Sterile pipettes supplied in peel-back packs printed with Lot No. and expiry date

## Product Applications

Fast nonvolumetric transfer of liquids

Drop counting

Blood banks

Hematology

Serology

Bacteriology

Immunology

Chemistry

General laboratory applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|              |                              |
|--------------|------------------------------|
| Art. No.     | E1414-0100                   |
| Pack Size    | 500 Pcs. (1 Pack × 500 Pcs.) |
| Length       | 160 mm                       |
| Diameter     | 5.0 mm                       |
| Sterile      | No                           |
| Autoclavable | No                           |
| Volume       | 1.0 ml                       |
| Tip Type     | standard                     |
| Material     | Polyethylene                 |
| Graduations  | Yes                          |
| Bulb draw    | 3.5 ml                       |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 1.0 ml Graduated Pasteur Pipette

A range of polyethylene pipettes which are manufactured to ISO standards. Ideal for aliquots, pasteur pipettes transfer small amounts of liquid with precision

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)