




1.0ml Micro Fine Tip Pasteur Pipette (Sterile), Ind. Wrapped



Product Highlights

- Low affinity surface
- Attached bulb avoids cross-contamination
- One-step biohazard disposal
- Shatter-proof low-density polyethylene, which is nontoxic and inert
- Can be frozen in liquid nitrogen
- Sterile pipettes supplied in peel-back packs printed with Lot No. and expiry date

Product Applications

Fast nonvolumetric transfer of liquids

Drop counting

Blood banks

Hematology

Serology

Bacteriology

Immunology

Chemistry

General laboratory applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







General Data

Art. No.	E1414-1111
Pack Size	400 Pcs. (1 Pack × 400 Pcs.)
Length	104 mm
Diameter	5.9 mm
Sterile	Yes
Autoclavable	No
Volume	1.0 ml
Тір Туре	Fine Tip
Material	Polyethylene
Graduations	No
Bulb draw	1.0 ml

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







More informations about 1.0ml Micro Fine Tip Pasteur Pipette (Sterile), Ind. Wrapped

A range of polyethylene pipettes which are manufactured to ISO standards. Ideal for aliquots, pasteur pipettes transfer small amounts of liquid with precision

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD 5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

