



7.0 ml Polystyrene Bijou Container (Sterile), without Label



Product Highlights

- Ultra clear sterile polystyrene container
- Centrifugation rate: 7.200 x g
- Temperature use: -10 °C to 70 °C
- 50 mm high x 23 mm diam (incl. cap)

Product Applications

Centrifugation

Sample preparation

Preparation of master mixes and buffers

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	E1412-0710
Overall Dimensions	50 mm high x 23 mm diameter
Sterile	Yes
Autoclavable	No
Volume	7.0 ml
Pack Size	700 Pcs. (1 Pack x 700 Pcs.)
Label Type	No Label
Material	Polystyrene
Max. temp.	70 °C
Temperature range	-10 to 70 °C
Temperature min. (SL)	-10 °C
g-Force max.	7.200 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 7.0 ml Polystyrene Bijou Container (Sterile), without Label





Ultra clear, sterile polystyrene bijou container with white polypropylene cap. Centrifugation rate 7200 x g. Temperature use: -10 °C to 70 °C.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 7.0 ml Polystyrene Bijou Container (Sterile) with Label Label Type: Plain Label  Create PDF Data Sheet	700 Pcs. (1 Pack × 700 Pcs.)	E1412-0711
 7.0 ml Polystyrene Bijou Container (Sterile) without Label Label Type: No Label  Create PDF Data Sheet	700 Pcs. (1 Pack × 700 Pcs.)	E1412-0710

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk