



## 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (sterile)



### Product Highlights

- **Graduated** at 10 µl, 50 µl and 100 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified** RNase, DNase, DNA and Pyrogen free.

## Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	S1113-1216
Pack Size	960 Tips (1 Stack rack × 960 Tips)
Color Code Wafer	Yellow
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Tip Type	UltraPoint®
Material	Polypropylene
Graduations	Yes
Color	Yellow
For use with	<p><b>Abimed:</b> BioPette-E, Discovery Comfort, LABMATE Optima</p> <p><b>Biohit®:</b> eLine, ePet, mLine, Picus, Proline</p> <p><b>BRAND®:</b> Transferpette, Transferpette Electronic, Transferpette S</p> <p><b>Capp®:</b> CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel</p> <p><b>Eppendorf®:</b> Reference, Reference 2, Research, Research plus, Research pro, Xplorer</p> <p><b>Gilson®:</b> Pipetman Concept, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p><b>Hamilton®:</b> SoftGrip Pipette, SofTouch</p> <p><b>HTL:</b> Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p><b>Nichiryo®:</b> 00-NP7, 50F, 50V, NLE, NPLO, NPX</p> <p><b>Rainin®:</b> Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p><b>Sartorius®:</b> eLine, ePet, mLine, Picus, Picus Nxt, Proline, Proline Plus, Tacta</p> <p><b>Socorex®:</b> Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital</p> <p><b>STARLAB®:</b> ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p><b>Thermo Fisher Scientific Inc.:</b> Finn Focus, Finn Focus Micro, Finn Focus Univ., Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Turquoise, Finnpipette Novus Yellow</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (sterile)













### UltraPoint orifice

Thin-wall orifice helps prevent sample hang up and makes it easier to touch off the last drop of sample

**The compact Stack Rack** provides 960 tips in one convenient unit, making an excellent choice for high throughput or multi-channel use. The wide base provides stability during use and the transparent cover protects the tips while not in use. The Stack Rack is autoclavable and all components are made from polypropylene and can be recycled.















## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1113-1806
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1113-1806-C
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1113-1816
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1113-1816-C
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1113-1706
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 9,600 Tips (100 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (100 Refills × 96 Tips)	S1113-1706-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1113-1716
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Refills × 96 Tips)	S1113-1716-C
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (1 Stack rack × 960 Tips)	S1113-1206
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack racks (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (10 Stack racks × 960 Tips)	S1113-1206-C
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (1 Stack rack × 960 Tips)	S1113-1216
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack racks (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (10 Stack racks × 960 Tips)	S1113-1216-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Bag (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	S1113-1006
 <a href="#">200 µl UltraPoint® graduated TipOne® Tip, Yellow, Bags (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 25,000 Tips (25 Bags × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	25,000 Tips (25 Bags × 1000 Tips)	S1113-1006-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)