



30 ml Polystyrene Universal Container, Plain Label



Product Highlights

- Conical base for concentration of particulate matter
- Aseptically produced in a fully automated and bio-controlled environment.
- Plain label
- Dimensions: 93 mm high x 31 mm diameter (including cap)
- Leak tested to 97kPa
- White HDPE caps with full width liner seal

Product Applications

Sample preparation

Centrifugation

Sample storage

Preparation of master mixes and buffers

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E1412-3011
Pack Size	400 Tubes (1 Bag x 400 Tubes)
Overall Dimensions	93 mm high x 31 mm diameter
Sterile	No
Autoclavable	No
Volume	30 ml
Label Type	Plain Label
Material	Polystyrene
Max. temp.	70 °C
Temperature range	-10 to 70 °C
Min. temp.	-10 °C
g-Force max.	1.600 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk















More informations about 30 ml Polystyrene Universal Container, Plain Label

Aseptically* produced crystal-clear polystyrene universal containers with leak proof screw cap. Centrifugation rate: 3,000 x g. Temperature use: -10 °C to 70 °C.

*Aseptically produced in a fully automated and bio-controlled environment.



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 30 ml Polystyrene Universal Container, No Label Material: Polystyrene Label Type: No Label  Create PDF Data Sheet	400 Tubes (1 Bag × 400 Tubes)	E1412-3010
 30 ml Polystyrene Universal Container, Printed Label Material: Polystyrene Label Type: Printed Label  Create PDF Data Sheet	400 Tubes (1 Bag × 400 Tubes)	E1412-3012
 30 ml Polypropylene Universal Container, No Label Material: Polypropylene Label Type: No Label  Create PDF Data Sheet	400 Tubes (1 Bag × 400 Tubes)	E1412-3020
 30 ml Polypropylene Universal Container, Plain Label Material: Polypropylene Label Type: Plain Label  Create PDF Data Sheet	400 Tubes (1 Bag × 400 Tubes)	E1412-3021
 30 ml Polystyrene Universal Container, Plain Label Material: Polystyrene Label Type: Plain Label  Create PDF Data Sheet	400 Tubes (1 Bag × 400 Tubes)	E1412-3011
 30 ml Polypropylene Universal Container, Printed Label Material: Polypropylene Label Type: Printed Label  Create PDF Data Sheet	400 Tubes (1 Bag × 400 Tubes)	E1412-3022

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk