



60 ml Polystyrene Container, Metal Cap, Printed Label (sterile)



Product Highlights

- Crystal clear polystyrene container with printed patient label
- 60ml
- Sterile
- Aluminium screw cap
- Temperature use: -10 °C to 70 °C.
- Dimensions: 60 mm high x 45 mm diameter (including cap)

Product Applications

Sample preparation

Sample storage

Preparation of master mixes and buffers

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E1412-6000
Pack Size	300 Pcs. (1 Box × 300 Pcs.)
Diameter	45 mm
Sterile	Yes
Autoclavable	No
Volume	60 ml
Label Type	Printed Label
Material	Polystyrene
Cap material	Aluminium
Max. temp.	70 °C
Temperature range	-10 to 70 °C
Min. temp.	-10 °C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 60 ml Polystyrene Container, Metal Cap, Printed Label (sterile)









Sterile, crystal-clear polystyrene container with metal cap and printed patient label. Temperature use: -10 °C to 70 °C.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 60 ml Polystyrene Container, Metal Cap, Printed Label (sterile) Volume: 60 ml  Create PDF Data Sheet	300 Pcs. (1 Box × 300 Pcs.)	E1412-6000
 150 ml Polystyrene Container, Metal Cap, non-sterile Volume: 150 ml  Create PDF Data Sheet	120 Pcs. (1 Box × 120 Pcs.)	E1412-1500
 100 ml Polystyrene Container, Metal Cap, non-sterile Volume: 100 ml  Create PDF Data Sheet	160 Pcs. (1 Box × 160 Pcs.)	E1412-1000
 250 ml Polystyrene Container, Metal Cap, non-sterile Volume: 250 ml  Create PDF Data Sheet	48 Pcs. (1 Box × 48 Pcs.)	E1412-2500

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk