



0.5 ml Homopolymer Microcentrifuge Tube, Natural



Product Highlights

- Ultra clear polypropylene tubes
- Firm and positive cap closure
- Suitable for boiling
- Rated at over 20,000 xg
- Certified RNase, DNase, DNA and PCR Inhibitor free

Product Applications

Sample preparation

Centrifugation

Sample storage

Sample boiling or freezing

Certified RNase, DNase, DNA and PCR Inhibitor free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I1405-1500
Pack Size	1,000 Tubes (1 Bag x 1000 Tubes)
Autoclavable	Yes
Volume	0.5 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 0.5 ml Homopolymer Microcentrifuge Tube, Natural

Polypropylene tubes with smaller lid footprint for better fit inside sample storage boxes. Amber tubes for light sensitive samples. Rated to over 20,000 x g. Temp. range: -80 °C to 121 °C. Autoclavable. Certified RNase, DNase, DNA and PCR Inhibitor free.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	0.5 ml Homopolymer Microcentrifuge Tube, Amber Color: Amber  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1507
	0.5 ml Homopolymer Microcentrifuge Tube, Natural Color: Natural  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1500
	0.5 ml Homopolymer Microcentrifuge Tube, Mixed Color: Mixed  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1508
	0.5 ml Homopolymer Microcentrifuge Tube, Yellow Color: Yellow  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1506
	0.5 ml Homopolymer Microcentrifuge Tube, Red Color: Red  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1504
	0.5 ml Homopolymer Microcentrifuge Tube, Green Color: Green  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1502
	0.5 ml Homopolymer Microcentrifuge Tube, Blue Color: Blue  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	I1405-1501

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk