



## 1.5 ml Homopolymer Microcentrifuge Tube, no Cap, Natural



### Product Highlights

- Ultra clear polypropylene tubes
- Suitable for boiling
- Rated to over 20,000 xg
- Certified RNase, DNase, DNA and PCR Inhibitor free

## Product Applications

Sample preparation

Centrifugation

Sample storage

Sample boiling or freezing

Certified RNase, DNase, DNA and PCR Inhibitor free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|                    |                               |
|--------------------|-------------------------------|
| Art. No.           | E1415-2850                    |
| Autoclavable       | Yes                           |
| Pack Size          | 500 Tubes (1 Bag × 500 Tubes) |
| Material           | Polypropylene                 |
| DNA free           | Yes                           |
| DNase free         | Yes                           |
| PCR inhibitor free | Yes                           |
| RNase free         | Yes                           |
| Color              | Natural                       |
| g-Force max.       | 20,000 x g                    |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 1.5 ml Homopolymer Microcentrifuge Tube, no Cap, Natural

Capless tubes for convenient collection of fractions or mix reagents. Rated to over 20,000 x g. Temp. range: -80 °C to 121 °C. Autoclavable. Certified RNase, DNase, DNA and PCR Inhibitor free.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)