



## Aluminium StarSeal (PCR)



### Product Highlights

- Suitable for standard PCR
- Single-coated seal consisting of aluminium backing, coated on one side with a pressure sensitive acrylate adhesive.
- The seal is supplied on a brown paper release liner
- Suitable for short-term storage
- Compatible with aqueous solutions and organic solvents
- Easy to pierce with a needle, metal probe or plastic pipette tip
- For use with polypropylene, polystyrene and polycarbonate plates
- Certified RNase, DNase, DNA and Endotoxin free
- Suitable for use between  $-70^{\circ}\text{C}$  and  $110^{\circ}\text{C}$
- Overall dimensions: 141 mm x 79 mm. Working dimensions: 121 mm x 79 mm

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	E2796-9792
Overall Dimensions	141 mm x 79 mm
Working Dimensions	121 x 79 mm
Sterile	No
Pack Size	100 Pcs. (1 Box x 100 Pcs.)
Material	Aluminium
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Pyrogen free	No
RNase free	Yes
Color	Silver
Plate type compatibility	Polypropylene, Polystyrene, Polycarbonate
Max. temp.	110 °C
Temperature range	-70 to 110 °C
Temperature min. (SL)	-70 °C
Breathable	no
Easy to Pierce with Needle or Metal Probe	yes
Easy to Pierce with Pipette Tip	yes
Low 'Tack to Touch' Adhesive	no
Low Auto-Fluorescence	no
Opaque	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



Optically clear	no
Resistant to DMSO	yes
Seal	Self-adhesive
Short-Term Storage and Incubation	yes
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	no
Superior Optical Clarity	no

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Aluminium StarSeal (PCR)

Aluminium StarSeal Sealing Film with pressure sensitive acrylate adhesive. The seal is piercable and resistant to DMSO.

Can be used for standard PCR and short term storage.

Suitable for use with polypropylene, polystyrene and polycarbonate plates. The use of a seal applicator is advised; strength of sealing is dependant on application pressure.



Suitable for use between -70 °C and 110 °C.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Plate Seal Applicator</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2928-7355
 <a href="#">Microtitre Sealing Brayer</a>	1 Piece (1 Pack × 1 Piece)	E9127-2940

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)