



Polypropylene PCR Sealing Film Strips, Clear



Product Highlights

- Supplied on convenient, perforated sheets of four strips
- Removable positioning tabs on either end provide easy application and removal
- Universal fit and specially formulated adhesive provides an excellent seal for all plate types
- Temperature range: -40°C to +120°C
- Certified RNase, DNase and nucleic acid free

Product Applications

Plate sealing

Real-Time PCR / qPCR

Standard PCR

Short term storage & incubation

Long term storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2796-2850
Pack Size	200 Pcs. (1 Box × 200 Pcs.)
Overall Dimensions	97 x 19 mm
Working Dimensions	76 x 19 mm
Material	Polypropylene
DNA free	No
DNase free	Yes
Pyrogen free	No
RNase free	Yes
Color	Clear
Plate type compatibility	Polypropylene
Max. temp.	120 °C
Temperature range	-40 to 120 °C
Min. temp.	-40 °C
Easy to Pierce with Needle or Metal Probe	no
Easy to Pierce with Pipette Tip	no
Low 'Tack to Touch' Adhesive	no
Low Auto-Fluorescence	no
Opaque	no
Optically clear	yes
Resistant to DMSO	no
Seal	Self-adhesive

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Short-Term Storage and Incubation	yes
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	no



More informations about Polypropylene PCR Sealing Film Strips, Clear



Plates seals are quick to apply and to remove. Eight self-adhesive seals offered by STARLAB are suitable for standard PCR. The Xtra-Clear Advanced Polyolefin Seal is also suitable for Real-Time PCR and the seal's superior optical clarity yields consistent results. All seals are suitable for short-term storage and incubation.

Important note for Real-Time PCR (qPCR)

Please check that both the plate AND the sealing option are suitable for Real-Time PCR applications.



Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Plate Seal Applicator	5 Pcs. (1 Pack × 5 Pcs.)	I2928-7355
	Microtitre Sealing Brayer	1 Piece (1 Pack × 1 Piece)	E9127-2940

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk