



## 0.45 $\mu\text{m}$ Syringe Filter, Cellulose Acetate (Sterile), Orange, $\square$ 33 mm



### Product Highlights

- **Cellulose Acetate Membrane** for general filtration and binding trace proteins
- High burst pressure ensures safe operation
- Large filtration surface area allows you to filter your sample more easily
- Low hold-up volume provides maximum recovery of your solution
- Surfactant-free membrane for cleanliness in downstream applications
- Filter type printed on the filter for ease of identification
- Convenient, larger pack sizes
- Gamma sterilised (SAL  $10^{-6}$ )

## Product Applications

0.22  $\mu\text{m}$  membranes for sterilization, sterility testing, “final filtration”

0.45  $\mu\text{m}$  membranes for analysis of bacteria, sterility testing, clarification and prefiltration

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E4780-1453
Pack Size	100 Pcs. (1 Bag × 100 Pcs.)
Diameter	33 mm
Sterile	Yes
Color	Orange
Membrane	Cellulose Acetate Membrane
Pore size	0.45 µm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)







## More informations about 0.45 $\mu\text{m}$ Syringe Filter, Cellulose Acetate (Sterile), Orange, $\varnothing$ 33 mm

### Nearly half the hold up volume of some other syringes on the market!

STARLAB's Syringe Filters are an economical choice due to the maximum recovery of your solution. Manufactured using the latest technology, these 33 $\varnothing$ mm diameter filters are suitable for volumes 10–200 $\varnothing$  ml. Use for cell culture media and additives, biological solutions and buffers.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.22 µm Syringe Filter, Cellulose Acetate (Sterile), Orange, Ø 33 mm</a> <b>Pore size:</b> 0.22 µm  <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Bag × 100 Pcs.)	E4780-1223
 <a href="#">0.45 µm Syringe Filter, Cellulose Acetate (Sterile), Orange, Ø 33 mm</a> <b>Pore size:</b> 0.45 µm  <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Bag × 100 Pcs.)	E4780-1453

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)