



200 µl XL Tip, Racked (non-sterile)



Product Highlights

- Compatible with a wide range of pipettes
- Made from polypropylene
- For some tips sterile racks or bags available
- Compatible with a wide range of pipettes
- Racked
- For some tips sterile racks or bags available

Product Applications

Pipetting liquids from and into long tubes
Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I1011-0920
Pack Size	1,632 Tips (8 Racks × 204 Tips)
Sterile	No
Filter tip	Non-filter
Volume range	1 - 200 µl
Volume	200 µl
Material	Polypropylene
DNase free	Yes
RNase free	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 200 µl XL Tip, Racked (non-sterile)











Extra long tips for easy handling of liquids in long tubes. Phase extraction and removal of lower phases.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl XL Tip Racked (sterile) Volume: 200 µl Sterile: Yes Packaging Unit: Racks  Create PDF Data Sheet	1,632 Tips (8 Racks × 204 Tips)	E1011-0921
 1,000 µl XL Bevelled Tip, Loose (non-sterile) Volume: 1000 µl Sterile: No Packaging Unit: Bag  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	E1011-1002
 200 µl XL Tip, Loose (non-sterile) Volume: 200 µl Sterile: No Packaging Unit: Bag  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	I1011-0020
 200 µl XL Tip, Racked (non-sterile) Volume: 200 µl Sterile: No Packaging Unit: Racks  Create PDF Data Sheet	1,632 Tips (8 Racks × 204 Tips)	I1011-0920
 1,000 µl XL Bevelled Tip, Racked (non-sterile) Volume: 1000 µl Sterile: No Packaging Unit: Racks  Create PDF Data Sheet	1,000 Tips (10 Racks × 100 Tips)	E1011-1430

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk