



50 ml Centrifuge Tube, Skirted, Loose (Sterile)



Product Highlights

- Clear, sterile polypropylene tubes with green HDPE screw caps
- White writing surface and clear, printed graduations
- Perfect for centrifugation of cells: USP Class VI, non-cytotoxic and sterile ($SAL 10^{-6}$) material
- Certified RNase, DNase, DNA, heavy-metal and Pyrogen free
- Max. centrifugation rate: 6,000 x g

Product Applications

Centrifugation

Sample preparation

Preparation of master mixes and buffers

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E1450-0400
Pack Size	500 Tubes (20 Bags × 25 Tubes)
Dimensions (W × H)	34.9 mm × 115.5 mm
External width	34.9 mm
Sterile	Yes
Autoclavable	Yes
Volume	50 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Clear
g-Force max.	12,000 x g
Non-cytotoxic	Yes
Tube Base	Skirted

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk









More informations about 50 ml Centrifuge Tube, Skirted, Loose (Sterile)

STARLAB's new centrifuge tubes are designed for the centrifugation of cells, but can also be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. The sterility of assurance level (SAL) of 10^{-6} fulfills highest demands on sterility, meeting EN 556-1 standard (Sterilisation of medical devices). Racked tubes are delivered in recyclable cardboard racks.

Maximum centrifugation rates: Conical Tubes 15 & 50 ml: 12,000 x g and Skirted Tube 50 ml: 6,000 x g



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	50 ml Centrifuge Tube, Conical, Loose (Sterile) Tube Base: Conical Packaging Unit: Bags  Create PDF Data Sheet	500 Tubes (20 Bags × 25 Tubes)	E1450-0100
	50 ml Centrifuge Tube, Conical, Racked (Sterile) Tube Base: Conical Packaging Unit: Racks  Create PDF Data Sheet	500 Tubes (20 Racks × 25 Tubes)	E1450-0700
	50 ml Centrifuge Tube, Skirted, Loose (Sterile) Tube Base: Skirted Packaging Unit: Bags  Create PDF Data Sheet	500 Tubes (20 Bags × 25 Tubes)	E1450-0400

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk