




StarTags Laboratory Labels (9.5 mm Diameter spot), orange



Product Highlights

- Suitable for autoclaving, boiling and freezing down to -85 °C
- Sticks permanently to all forms of tubes, plates, vials and bottles
- No sticky residue left following removal

Product Applications

Labelling

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E9185-0003
Pack Size	2,000 Pcs. (1 Box × 2000 Pcs.)
Diameter	9.5 mm
Autoclavable	Yes
Color	Orange
Operating temperature	Suitable for boiling and freezing down to -85 °C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about StarTags Laboratory Labels (9.5 mm Diameter spot), orange













StarTags: One label you'll want to stick with!

Pre-cut labels that will accept any form of marker pen and will stick permanently in most conditions. StarTags are suitable for use in the freezer, centrifuge, waterbath, thermal cycler, shaker or vortexer. Laser sheets are available in some options, allowing you print out large numbers of labels which is both time saving and convenient.

NOTE: Labels must be applied to dry, unfrozen samples at room temperature. StarTags are not suitable for use with low-retention/siliconised tubes. A4 Laser Sheets are supplied with set-up instructions for MS Word and can also be downloaded from the STARLAB website. Sheet labels are only suitable for use with laser printers and can not be used with bubble jet printers.



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	StarTags Laboratory Labels (9.5 mm Diameter spot), blue Color: Blue	2,000 Pcs. (1 Box × 2000 Pcs.)	E9185-0001
	 Create PDF Data Sheet		
	StarTags Laboratory Labels (9.5 mm Diameter spot), green Color: Green	2,000 Pcs. (1 Box × 2000 Pcs.)	E9185-0002
	 Create PDF Data Sheet		
	StarTags Laboratory Labels (9.5 mm Diameter spot), orange Color: Orange	2,000 Pcs. (1 Box × 2000 Pcs.)	E9185-0003
	 Create PDF Data Sheet		
	StarTags Laboratory Labels (9.5 mm Diameter spot), red Color: Red	2,000 Pcs. (1 Box × 2000 Pcs.)	E9185-0004
	 Create PDF Data Sheet		
	StarTags Laboratory Labels (9.5 mm Diameter spot), yellow Color: Yellow	2,000 Pcs. (1 Box × 2000 Pcs.)	E9185-0006
	 Create PDF Data Sheet		
	StarTags Laboratory Labels (9.5 mm Diameter spot), white Color: White	2,000 Pcs. (1 Box × 2000 Pcs.)	E9185-0009
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk