



## Starlab ID Tape Core Inserts, Pack of 5



### Product Highlights

- Core inserts to hold the 55 m long tape in the dispenser available separately
- Core inserts for 55 m long tape available separately

## Product Applications

Labelling

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E9026-2005
Pack Size	5 Pcs. (1 Box × 5 Pcs.)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Starlab ID Tape Core Inserts, Pack of 5

### Dispense with disorganisation!

- › Custom-made acrylic dispenser which holds both lengths of tape
- › Suction feet ensures the dispenser stays firmly on the bench, even when using one hand to dispense the tape
- › Core inserts to hold the 55 m long tape in the dispenser available separately

Dispenser dimensions (WxDxH): 205 x 185 x 125 mm