



## 1L Distel DS Disinfectant Concentrate



### Product Highlights

- 1 L Bottle Distel DS Concentrate (includes 1 x 5.0 ml Dosing Head)
- The new-style disinfectant system to replace wipes
- Effective against bacteria, yeast and viruses in just 3 minutes
- 1 L of concentrate provides 2000 x 5.0 ml doses
- Volume markers on bottles aid quick and easy dilution
- Dosing heads provide consistent dosage, every time
- Save plastic waste! Reuse the bottles again and again

## Product Applications

Disinfection  
Decontamination  
Cleaning  
Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|           |                              |
|-----------|------------------------------|
| Art. No.  | CCH050201                    |
| Volume    | 1 l                          |
| Pack Size | 1 Piece (1 Bottle × 1 Piece) |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)





## More informations about 1L Distel DS Disinfectant Concentrate

**The new-style disinfectant system to replace wipes!** Take advantage of this new system where you can use any cloth/towel and dose it with the amount you need. Distel DS is a convenient and easy to use surface disinfectant for the cleaning and disinfection of hard surfaces and non-invasive medical devices. Simply dilute 1 litre of concentrate with tap water to make 10 litres of disinfectant that is effective against bacteria, yeast and viruses in just 3 minutes! The clever dosing head on the bottles ensures a consistent 5.0 ml dose of disinfectant every time. A starter pack is available that contains all you need to get going (1L disinfectant concentrate and five reusable dosing bottles); you can then purchase the 1L concentrate on its own so you can reuse the bottles over and again. They're made from polypropylen so they're recyclable too. The whole system is easy, green and results in less plastic waste. Eliminate single-use plastic wipes now with Distel DS - just dilute, dose and wipe using a towel or cloth of your choice.



## All Variations

| PRODUCT NAME   | PACKAGING SIZE               | ART. NO.  |
|--|------------------------------|-----------|
| <br><a href="#">1L Distel DS Disinfectant Concentrate</a><br><b>Packaging Unit:</b> Bottle<br> <a href="#">Create PDF Data Sheet</a> | 1 Piece (1 Bottle × 1 Piece) | CCH050201 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)