



Autoclave bags 415 x 600 mm



Product Highlights

- Polypropylene bags for autoclaving up to 136°C
- 48 microns
- Printed with the Biohazard Symbol
- Use for contaminated waste disposal
- White write-on panel

Product Applications

Disposal

Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	E3300-4070
Pack Size	200 Pcs. (1 Box × 200 Pcs.)
Dimensions (W × H)	415 mm × 600 mm
External width	415 mm
Autoclavable	Yes
Material	Polypropylene
Color	Clear

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Autoclave bags 415 x 600 mm







Autoclave bags

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Autoclave bags 310 x 660 mm Dimensions (W x H): 310 mm x 660 mm	200 Pcs. (1 Box x 200 Pcs.)	E3300-3166
	Create PDF Data Sheet		
	Autoclave bags 415 x 600 mm Dimensions (W x H): 415 mm x 600 mm	200 Pcs. (1 Box x 200 Pcs.)	E3300-4070
	Create PDF Data Sheet		
	Autoclave bags 610 x 810 mm Dimensions (W x H): 610 mm x 810 mm	200 Pcs. (1 Box x 200 Pcs.)	E3300-6078
	Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk