



200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile), Case



Product Highlights

- **Graduated** at 10 µl, 50 µl and 100 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified** RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids
Dispensing liquids
Mixing liquids
Filling plates and reaction vessels
Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	S1113-1716-C
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1 - 200 µl
Volume	200 µl
Pack Size	7,680 Tips (80 Refills × 96 Tips)
Color Code Wafer	Yellow
Tip Type	UltraPoint®
Material	Polypropylene
Graduations	Yes
Color	Yellow
For use with brand + product (SL)	<p>Abimed: BioPette-E, Discovery Comfort, LABMATE Optima</p> <p>Biohit®: eLine, ePet, mLine, Picus, Proline</p> <p>BRAND®: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S</p> <p>Capp®: CappAero Comfort, CappAero Multi, CappAero Universal, Cappronic Multi, Cappronic Single, ecopipette, Singlechannel</p> <p>Eppendorf®: Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer</p> <p>Gilson®: Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p>Hamilton®: SoftGrip Pipette, SofTouch</p> <p>HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p>Nichiryo®: 00-NP7, 50F, 50V, NLE, NPLO, NPX</p> <p>Rainin®: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p>Sartorius®: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta</p> <p>Socorex®: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital</p> <p>STARLAB®: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p>Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finn Focus Univ., Finn pipette Digital, Finn pipette F1/F2, Finn pipette Novus, Finn pipette Novus Turquoise, Finn pipette Novus Yellow</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile), Case

Experience effortless precision with [TipOne® Pipette Tips](#)

Elevate your daily lab work with [TipOne® pipette tips](#)! Enjoy precise, reliable pipetting every time—no more wasted samples or frustrating fit issues. [TipOne®](#) tips are compatible with all popular pipettes for accurate results you can trust, every single day.

Why [TipOne®](#) makes labs better:

Universal tip fit: Simplify your bench and inventory. The tip collar of [TipOne®](#) tips is engineered to fit all popular pipettes without any loss of accuracy or precision.

Made for precision: Expert design and the use of high-quality virgin materials means our tips provide consistent accuracy to help you get reliable results and reduce costly rework.

Pure quality: Sterile [TipOne®](#) tips are validated according to EN ISO 11137-2:2007. All [TipOne®](#) tips are certified RNase, DNase, DNA and pyrogen free.

Flexible Packaging:

Racked: Sturdy, reusable racks stay stable on your bench. Quickly reload with [TipOne®](#) refills for reduced plastic waste and cost. Sterile options available.

Refill: Load tips into your racks with a few easy steps. Refills are available sterile or non-sterile, helping you cut down on plastic waste (and expenses).

Stack Rack: Handle large experiments effortlessly—dispense 960 tips from one compact unit that is perfect for multi-channel pipettes and busy workflows. Sterile options available.

Loose: For teaching or routine work, stock up on affordable quality tips. Supplied in convenient resealable bags.

Join our journey to sustainable pipetting!

Reduce your lab's plastic waste! [TipOne®](#) racks and refill system components use recycled materials and are fully recyclable (except shrink wrap, loose bag, and tips) via the [TipOne®](#) recycling service.













We'll collect your used [TipOne®](#) plastic, recycle it and use the recycled plastic to make more [TipOne®](#) products. No more guessing if or how your recycling gets reused—[TipOne®](#) closes the loop for you!

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk















All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (10 Racks × 96 Tips)	S1113-1806
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	7,680 Tips (80 Racks × 96 Tips)	S1113-1806-C
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (10 Racks × 96 Tips)	S1113-1816
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	7,680 Tips (80 Racks × 96 Tips)	S1113-1816-C
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (10 Refills × 96 Tips)	S1113-1706
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p>	9,600 Tips (100 Refills × 96 Tips)	S1113-1706-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1113-1716
 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Refills × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Refills × 96 Tips)	S1113-1716-C
 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)  Create PDF Data Sheet	960 Tips (1 Stack rack × 960 Tips)	S1113-1206
 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)  Create PDF Data Sheet	9,600 Tips (10 Stack racks × 960 Tips)	S1113-1206-C
 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (sterile) Sterile: Yes Pack Size: 960 Tips (1 Stack rack × 960 Tips)  Create PDF Data Sheet	960 Tips (1 Stack rack × 960 Tips)	S1113-1216
 200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack racks (sterile), Case Sterile: Yes Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)  Create PDF Data Sheet	9,600 Tips (10 Stack racks × 960 Tips)	S1113-1216-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)</p> <p> Create PDF Data Sheet</p>	1,000 Tips (1 Bag × 1000 Tips)	S1113-1006
 <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Bags (non-sterile), Case Sterile: No Pack Size: 25,000 Tips (25 Bags × 1000 Tips)</p> <p> Create PDF Data Sheet</p>	25,000 Tips (25 Bags × 1000 Tips)	S1113-1006-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk