



## 20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Racks (sterile), Case



### Product Highlights

- Graduated at 10 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal tip fit.
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- Racks can be retained and reused with TipOne® refills
- Rack base made with ~25 % recycled plastic
- All components, excluding tips and shrink wrap, are recyclable
- Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

## Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	S1123-1810-C
Pack Size	7,680 Tips (80 Racks × 96 Tips)
Color Code Wafer	Yellow
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	1.0 - 20 µl
Volume	20 µl
Tip Type	UltraPoint®
Material	Polypropylene
Graduations	Yes
Color	Natural
Sustainability attributes	<p>Racks can be retained and reused with TipOne® refills.</p> <p>Rack base made with ~25 % recycled plastic.</p> <p>All components, excluding film wrap and tips, are recyclable.</p> <p>Buy by the case to save on packaging and transportation.</p> <p>Reuse the case box to send back your recycling.</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



For use with

**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima

**Biohit:** eLine, mLine, Picus, Proline

**BRAND:** Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

**Capp:** CappAero Comfort, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel

**Eppendorf:** Reference, Reference 2, Research, Research plus

**Gilson:** Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

**HTL:** Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

**Nichiryo:** 50F, NLE, NPLO, NPX

**Rainin:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

**Sartorius:** eLine, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

**Socorex:** Acura 825, Calibra digital

**STARLAB:** ErgoOne, ErgoOne E, MicroOne, StarPet

**Thermo Fisher Scientific Inc.:** Finn Focus, Finn Focus Univ., Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Yellow

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Racks (sterile), Case

### Nothing goes through!

Protect your samples from unwanted contamination with TipOne® filter tips. Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your pipette, but they can contaminate other samples later. TipOne® filter tips block over 99% of aerosols\*, making them one of the best filters available.

Sterile TipOne® tips are validated according to EN ISO 11137-2:2007.

### Join our journey to sustainable pipetting.









TipOne® tip racks are built to last. When you have used your filter tips, the rack can be reused again and again with the cost-effective TipOne® refills (non-sterile, sterile and sterile filter tip options available). TipOne® racks have been tested and they can withstand at least 100 autoclave cycles (rack with wafer inside) without damage.

The base of the rack is made with 25 % recycled plastic. All components, excluding tips and shrink, are recyclable. Sign up for our TipOne® recycling service!

\*See TipOne® brochure in downloads section for details.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Refills (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1123-1710
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Refills (sterile), Case</a> <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Refills × 96 Tips)	S1123-1710-C
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Racks (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1123-1810
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Racks (sterile), Case</a> <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1123-1810-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)