



## 1,000 µl TipOne® Filter Tip, Natural, Racks (sterile), Case



### Product Highlights

- 100 µl, 250 µl, 500 µl and 1000 µl
- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- Certified RNase, DNase, DNA and Pyrogen free

## Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	S1126-7810-C
Color Code Wafer	Blue
Material	Polypropylene
Pack Size	3,840 Tips (40 Racks × 96 Tips)
Autoclavable	No
Sterile	Yes
Graduations	No
Filter tip	Filter
Volume range	100 - 1,000 µl
Color	Natural
Volume	1000 µl
For use with	<b>Abimed</b> : BioPette-E, Discovery Comfort, LABMATE Optima <b>Biohit</b> : ePet, mLine, Proline <b>BRAND</b> : Transferpette, Transferpette Electronic, Transferpette S <b>Capp</b> : CappAero Comfort, CappAero Universal, Capptronic Single, ecopipette, Singlechannel <b>Eppendorf</b> : Reference, Reference 2, Research, Research plus, Research pro, Xplorer <b>Gilson</b> : Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra <b>Hamilton</b> : SoftGrip Pipette, SofTouch <b>HTL</b> : Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette <b>Nichiryo</b> : 31N, 50F, 50V, NLE, NPLO, NPX <b>Rainin</b> : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus <b>Sartorius</b> : ePet, mLine, Proline, Proline Plus, Tacta <b>Socorex</b> : Acura 815 manual, Acura 825, Acura electro 925 <b>STARLAB</b> : ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E <b>Thermo Fisher Scientific Inc.</b> : Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 1,000 µl TipOne® Filter Tip, Natural, Racks (sterile), Case

### Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.









During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.

TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

*\*The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1,000 µl TipOne® Filter Tip, Natural, Racks (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1126-7810
 <a href="#">1,000 µl TipOne® Filter Tip, Natural, Racks (sterile), Case</a> <b>Pack Size:</b> 3,840 Tips (40 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	3,840 Tips (40 Racks × 96 Tips)	S1126-7810-C
 <a href="#">1,000 µl TipOne® Filter Tip, Natural, Refills (sterile), Case</a> <b>Pack Size:</b> 3,840 Tips (40 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	3,840 Tips (40 Refills × 96 Tips)	S1126-7710-C
 <a href="#">1,000 µl TipOne® Filter Tip, Natural, Refills (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1126-7710

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)