



1,000 µl TipOne® Filter Tip, Natural, Racks (sterile), Case



Product Highlights

- Highly effective HDPE filter that contains no additives to interfere with samples
- Universal tip fit
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- Racks can be retained and reused with TipOne refills
- All components, excluding tips and shrink wrap, are recyclable
- Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

Product Applications

DNA applications (e.g. PCR)
RNA applications (e.g. Gene expression analysis)
Protein applications (e.g. Antibody Research)
Cell Culture applications (e.g. Media)
Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD







General Data

Art. No.	S1126-7810-C		
Pack Size	3,840 Tips (40 Racks × 96 Tips)		
Color Code Wafer	Blue		
Sterile	Yes		
Autoclavable	No		
Filter tip	Filter		
Volume range	100 - 1,000 μl		
Volume	1000 μl		
Тір Туре	standard		
Material	Polypropylene		
Graduations	No		
Color	Natural		
Sustainability attributes	Racks can be retained and reused with TipOne [®] refills. Rack base made with ~25 % recycled plastic. All components, excluding tips and shrink wrap, are recyclable. Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.









For use with

Abimed: BioPette-E, Discovery Comfort, LABMATE Optima

Biohit*: ePet, mLine, Proline

BRAND°: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

Capp*: CappAero Comfort, CappAero Universal, Capptronic Single, ecopipette, Singlechannel

Eppendorf*: Reference, Reference 2, Research, Research plus, Research pro, Xplorer

Gilson°: Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

Hamilton^{*}: SoftGrip Pipette, SofTouch

HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labrate Pro, Optipette

Nichiryo°: 31N, 50F, 50V, NLE, NPLO, NPX

Rainin*: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus Sartorius*: ePet, mLine, Proline, Proline Plus, Tacta

Socorex^{*}: Acura 815 manual, Acura 825, Acura electro 925 STARLAB^{*}: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

Thermo Fisher Scientific Inc.: Finn Focus, Finnpipette Digital, Finnpipette F1/F2,

Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD







More informations about 1,000 µl TipOne® Filter Tip, Natural, Racks (sterile), Case

Nothing goes through!

Protect your samples from unwanted contamination with TipOne[®] filter tips. Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your pipette, but they can contaminate other samples later. TipOne[®] filter tips block over 99% of aerosols*, making them one of the best filters available.

Sterile TipOne® tips are validated according to EN ISO 11137-2:2007.

Join our journey to sustainable pipetting.

TipOne[®] tip racks are built to last. When you have used your filter tips, the rack can be reused again and again with the cost-effective TipOne[®] refills (non-sterile, sterile and sterile filter tip options available). TipOne[®] racks have been tested and they can withstand at least 100 autoclave cycles (rack with wafer inside) without damage.

The base of the rack is made with 25 % recycled plastic. All components, excluding tips and shrink, are recyclable. Sign up for our TipOne* recycling service!

*See TipOne® brochure in downloads section for details.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD 5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	1,000 µl TipOne® Filter Tip, Natural, Racks (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips) Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1126-7810
	1,000 µl TipOne® Filter Tip, Natural, Racks (sterile), Case Pack Size: 3,840 Tips (40 Racks × 96 Tips) Create PDF Data Sheet	3,840 Tips (40 Racks × 96 Tips)	S1126-7810-C
	1,000 µl TipOne® Filter Tip, Natural, Refills (sterile), Case Pack Size: 3,840 Tips (40 Refills × 96 Tips) Create PDF Data Sheet	3,840 Tips (40 Refills × 96 Tips)	S1126-7710-C
=	1,000 µl TipOne® Filter Tip, Natural, Refills (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips) Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1126-7710

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.



