



10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile), Case



Product Highlights

- Graduated at 1.0 µl, 5.0 µl, 10 µl and 20 µl
- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- Certified RNase, DNase, DNA and Pyrogen free.

Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	S1120-3710-C
Pack Size	7,680 Tips (80 Refills × 96 Tips)
Color Code Wafer	Red
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	0.1 - 20 µl
Volume	20 µl
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with	<p>Abimed: BioPette-E, Discovery Comfort, LABMATE Optima</p> <p>Biohit®: eLine, ePet, mLine, Picus, Proline</p> <p>BRAND®: Transferpette Electronic, Transferpette S</p> <p>Capp®: CappAero Comfort, CappAero Multi, CappAero Universal, CappAero Universal short, ecopipette, Singlechannel</p> <p>Eppendorf®: Reference, Reference 2, Research, Research plus, Research pro, Xplorer</p> <p>Gilson®: Pipetman Concept, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p>Hamilton®: SoftGrip Pipette, SofTouch</p> <p>HTL: Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p>Nichiryo®: 00-NP7, 50F, NLE, NPLO, NPX</p> <p>Rainin®: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p>Sartorius®: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta</p> <p>Socorex®: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital</p> <p>STARLAB®: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p>Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finn pipette Digital, Finn pipette F1/F2, Finn pipette Novus, Finn pipette Novus Pink</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile), Case

Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.









During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.

TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

**The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips)	960 Tips (10 Refills × 96 Tips)	S1120-3710
	 Create PDF Data Sheet		
	10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile), Case Pack Size: 7,680 Tips (80 Refills × 96 Tips)	7,680 Tips (80 Refills × 96 Tips)	S1120-3710-C
	 Create PDF Data Sheet		
	10/20 µl XL graduated TipOne® Filter Tip, Racks (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips)	960 Tips (10 Racks × 96 Tips)	S1120-3810
	 Create PDF Data Sheet		
	10/20 µl XL graduated TipOne® Filter Tip, Racks (sterile), Case Pack Size: 7,680 Tips (80 Racks × 96 Tips)	7,680 Tips (80 Racks × 96 Tips)	S1120-3810-C
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk