



## Large Volume Tips 10 ml



### Product Highlights

- From 1200  $\mu$ l tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

## Product Applications

Pipetting of large volumes of liquids

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|              |                |
|--------------|----------------|
| Art. No.     | See variations |
| Autoclavable | Yes            |
| Filter tip   | Non-filter     |
| Volume       | 10 ml          |
| Material     | Polypropylene  |
| DNA free     | Yes            |
| DNase free   | Yes            |
| Pyrogen free | Yes            |
| RNase free   | Yes            |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Large Volume Tips 10 ml











[Large Volume Tips](#)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)





## All Variations

| PRODUCT NAME  | PACKAGING SIZE                | ART. NO.   |
|---|-------------------------------|------------|
|  <p><a href="#">10 ml Finn/Socorex-Style Blue Tip (non-sterile)</a><br/> <b>Design Tips:</b> Finn/Socorex-Style<br/> <b>Sterile:</b> No<br/> <b>Pack Size:</b> 100 Tips (1 Bag × 100 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>      | 100 Tips (1 Bag × 100 Tips)   | I1052-0000 |
|  <p><a href="#">10 ml Finn/Socorex-Style Blue Tip (non-sterile)</a><br/> <b>Design Tips:</b> Finn/Socorex-Style<br/> <b>Sterile:</b> No<br/> <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>    | 250 Tips (10 Racks × 25 Tips) | I1052-0001 |
|  <p><a href="#">10 ml Finn/Socorex-Style Blue Tip (sterile)</a><br/> <b>Design Tips:</b> Finn/Socorex-Style<br/> <b>Sterile:</b> Yes<br/> <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>   | 250 Tips (10 Racks × 25 Tips) | I1052-0010 |
|  <p><a href="#">10 ml Gilson-Style Graduated Macro Tip (non-sterile)</a><br/> <b>Design Tips:</b> Gilson-Style<br/> <b>Sterile:</b> No<br/> <b>Pack Size:</b> 100 Tips (1 Bag × 100 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>   | 100 Tips (1 Bag × 100 Tips)   | I1054-0000 |
|  <p><a href="#">10 ml Gilson-Style Graduated Macro Tip (non-sterile)</a><br/> <b>Design Tips:</b> Gilson-Style<br/> <b>Sterile:</b> No<br/> <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p> | 250 Tips (10 Racks × 25 Tips) | I1054-0001 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



| PRODUCT NAME  | PACKAGING SIZE                | ART. NO.   |
|---|-------------------------------|------------|
|  <p><a href="#">10 ml Gilson-Style Graduated Macro Tip (sterile)</a><br/><b>Design Tips:</b> Gilson-Style<br/><b>Sterile:</b> Yes<br/><b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p> | 250 Tips (10 Racks × 25 Tips) | I1054-0011 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)