




Flat Gel-Loading Tips



Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

Gel electrophoresis

Loading samples onto gels (eg. polyacrylamide or agarose)

Handling of small volumes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Sterile	No
Volume	200 µl
Pack Size	200 Tips (1 Rack × 200 Tips)
Tip Type	flat gel-loading

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Flat Gel-Loading Tips





All tips are certified RNase and DNase free. Packaging of these tips differs from the TipOne® range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl Flat Gel-Loading Tip (Tip Thickness 0.20 mm) Tip Thickness: 0.20 mm  Create PDF Data Sheet	200 Tips (1 Rack × 200 Tips)	I1022-2600
 200 µl Flat Gel-Loading Tip (Tip Thickness 0.40 mm) Tip Thickness: 0.40 mm  Create PDF Data Sheet	200 Tips (1 Rack × 200 Tips)	I1022-4600

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk